




**GT040-100
Disposable Surgical Blade**

Size: 10#,10A#,11#,12#,15#.20#,21#,22#,23#,24#
Material: Carbon steel or stainless steel
Packing: 1pc/Aluminium Foil Bag,
100or200pcs /box, 5000pcs/ctn




**GT042-100
Blood Lancet**

Stainless steel Blood lancet
Length: 40mm



**GT040-300
Skin Graft Blade**

Size: 158mm
Material : Carbon Steel Or Stainless Steel
Packing: 1pc/Aluminium Foil Bag,
100pcs /box, 5000pcs/ctn



**GT040-200
Stitich Cutter Blade**

Size: 65mm, 110mm
Material : Carbon Steel Or Stainless Steel
Packing:1pc/Aluminium Foil Bag,
100pcs /box, 5000pcs/ctn



**GT042-200
Plastic Handle Lancet--Twist Top**

Size	21G	23G	26G	28G	30G
Color	Orange	Green	Yellow	Pink	Purple


Sterilized by Gamma
Packing: 100pcs/box, 20000pcs/ctn

**GT042-202
Softclix Blood Lancet**

Plastic handle Flat Lancet

Size	28G	30G
Color	White	White

Sterilized by gamma
Packing:100pcs/box,20000pcs/ctn



**GT041-100
Disposable Scalpel**

Size:10#-36#
Material : Carbon Steel Or Stainless Steel
Packing:1pc/aluminium Foil Bag,
10pcs /box,500pcs/ctn



**GT041-200
Safety Surgical Scalpel**

Size:10#-36#
Material: Carbon Steel Or Stainless Steel,
Sterilized By Gamma
Packing:1pc/aluminium Foil Bag,
10pcs/box,500pcs/ctn



**GT043-302
Sterile Safety Lancet**

Size	21Gx2.2mm	23Gx2.3mm	26Gx1.8mm	28Gx1.8mm	30Gx1.8mm
Color	Orange	Green	Yellow	Pink	Purple



**GT042-400
Lancing Device**

Material : Plastic. 5 Depth Adjustment
Material : Plastic. 10 Depth Adjustment
Material : Plastic. 6 Depth Adjustment
Individually Packing: Box or Hard Blister Packaging.



**GT041-600
Surgical Handle**

Size: 3#/12.1cm, 4#/14cm
Material: Stainless Steel
Packing: 1pc/aluminium Foil Bag



Size	21G	26G	28G	30G
Color	Green	Yellow	Orange	Purple

**GT042-304
Safety Lancet**

Mini safety lancet, independent used type.
Length: 37.99mm
Width: 6.45mm
Depth: 1.8mm
Bleeding amount: 100-200µL



**GT043-301A
Pressure Activated Blood Lancet I**

Size	21Gx2.2mm	21Gx2.4mm	23Gx1.8mm	26Gx1.8mm	28Gx1.8mm	30Gx1.8mm
Color	Orange	Purple	Green	Yellow	Pink	Light Purple



**GT043-301B
Pressure Activated Blood Lancet II**



**GT042-303
Safety Lancet**

Size 18G 2.2mm blade 21G 2.4mm blade 21G 2.4mm 23G 2.2mm 26G 1.8mm 28G 1.8mm 30G 1.6mm
Color Pink Light Blue Green Blue Yellow Orange Purple

Size	18G 2.2mm Blade	21G 2.4mm Blade	21G 2.4mm	23G 2.2mm	26G 1.8mm	28G 1.8mm	30G 1.6mm
Color	Pink	Light Blue	Green	Blue	Yellow	Orange	Purple

Disposable Medical Razor

Carbon steel or stainless steel, single or double blade.
All razors come with the highest quality stainless steel blades or carbon steel, comfort-coated for smooth shaves.



GT040-200A



GT040-300



GT040-200B

GT040-200C



**GT043-100
Disposable Medical Razor**

Carbon steel or stainless steel, single or double blade.
All razors come with the highest quality Stainless Steel blades or carbon steel, Comfort-Coated for smooth shaves.
Packing: 100pcs/box, 1000pcs/carton, 63x32x21cm




Surgical Suture

Features:

- 1) Suture material: polyglycolic acid (synthetic absorbable sutures), chromic and plain catgut, braided silk, monofilament nylon, monofilament polypropylene, braided polyester, and monofilament stainless steel wire
 - 2) thread 45cm, 60cm, 75cm, 100cm, 125cm, 150cm ect.
 - 3) Thread diameter: 8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2, 3. ect
 - 4) Needle shape: 1/2 circle, 3/8 circle, 1/4 circle, 5/8 circle, straight cutting, J-shape
 - 5) Needle end type: rolled end, drilled end, regular-eye end
 - 4) Needle body: taper point, blunt taper point, cuttinh edge, reverse cutting edge, tapercut, micro-point spatula curved
 - 5) Sterilization method: gamma radiation or ethylene oxide (can be used directly without being re-sterilized before use)
- Inner packing: 1 set/sealed polyester and aluminum foil sachet

Outer packing:


12 foil sachets/printed paper box or plastic container 50boxes/carton(28X28X34cm)




Thread length:	45cm, 60cm, 75cm, 100cm, 125cm, 150cm etc.
Thread diameter:	8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2, 3 etc.
Needle length:	6mm, 8mm, 12mm, 18mm, 22mm, 30mm, 35mm, 40mm, 50mm etc.
Needle curvature:	straight, 1/2 circle, 1/2 circle(double), 1/4 circle, 1/4 circle(double), 3/8 circle, 3/8 circle(double), 5/8 circle, loop round
Cross section:	round bodied, round bodied(heavy), curved cutting, curved cutting(heavy), reverse cutting, reverse cutting(heavy), tapercut, micro-point spatula curved etc.




GT044-100
Chromic catgut




GT044-200
Plain catgut




GT044-700
Polyglycolic Acid




GT044-300
Silk braided



GT044-400
Nylon monofilament



GT044-500
Polyester braided



GT044-600
Polypropylene