

# Foley Catheter

Made from natural latex

Silicone coated

Available with the option of plastic valve or rubber valve

Colour-coded for visualisation of size

## GT029-610

### Closed Wound Drainage System (Spring)

1. This product includes 3-spring evacuator, PVC tubing, Y connector, PVC drainage tube and stainless steel trocar.
2. The design of 3-spring evacuator minimizes the Collapse-Resistance.
3. The PVC drainage tube inside can be replaced by our other silicone drains, if needed.

Size available : 200ml,400ml,800ml(8Fr, 10Fr,12Fr, 14Fr, 16Fr, 18Fr,20Fr,22Fr,24Fr)  
Packed :Blister Bag

## GT029-620

### Closed Wound Drainage System (Hollow)

1. Aseptic PE material.
  2. Hollow non-spring wound drainage reservoir with PVC drainage tube and trocar, bed sheet clamp and robot clamp.
- Size available : 200ml,400ml(8Fr, 10Fr,12Fr, 14Fr, 16Fr, 18Fr, 20Fr,22Fr,24Fr)  
Packing :Blister Bag



## GT069-112

### Silicone Reservoir

Made of medical grade silicone  
Provided with internal anti-reflux valve  
Bulb silicone walls provide easy activation of suction  
Regular system is to be used together with silicone flat drains, but we also recommend innovative solutions which are combined with our round fluted drains or round perforated drains.  
Size:100ml,150ml,200ml,400ml  
Packing :Blister bag




Round Fluted

Round Channel

## GT069-112T

### Silicone Round Channel/Fluted Drains

Made of medical grade silicone  
Utilize internal channels or fluted, to swiftly remove exudate Independent channels facilitate drainage and reduce risk of occlusion Made by one time, no connector ensures patients' comfort upon removal Radio-opaque line through the length for X-ray visualization  
Available with or without "Three Face" stainless steel trocar  
Size and length can be customized  
Size:FR 10,FR 15,FR 19,FR 24,FR 28  
Packing :Blister bag



Flat Perforated

Flat Fluted



## GT022-200 2-way Latex Foley Catheter

Made from natural latex

Silicone coated

Plastic valve or Rubber valve is available

Color-coated for visualization of size

Fr6-Fr10 (Pediatric)

Length:270mm Balloon Capacity:3cc/5cc

Fr12-Fr22 (Standard)

Length:400mm Balloon Capacity:5cc/10cc/30cc

Fr24-Fr26 (Standard)

Length:400mm Balloon Capacity:10cc/30cc

Packed in individual Blister bag

10pcs/box, 500pcs/ctn, Dimension:52X43X47.5 cm



## GT022-300

### 3-Way latex Foley Catheter

Made from natural latex

Silicone Coated

Plastic valve or Rubber valve is available

Used for clinical urine catheterization, injection and irrigation.

Length: 400mm

Size: Balloon Capacity

Fr16-Fr26 30cc

packing:Blister Bag



## GT022-400

### 2-Way Female Foley Catheter

Made from natural latex

Silicone Coated

Plastic valve or Rubber valve is available

Used for clinical urine catheterization of Female

Length: 260mm

Size: Balloon Capacity

Fr12-Fr26 30cc

packing:Blister Bag



## GT022-500

### 2-Way Tiemann Foley Catheter

Made from natural latex

Silicone Coated

Plastic valve or Rubber valve is available

Used for clinical urine catheterization of male


Length: 260mm

Size: Balloon Capacity

Fr12-Fr26 30cc

packing:Blister Bag





**GT022-100  
Urethral Catheter**

Made from natural latex  
Silicone Coated  
Plastic valve or Rubber valve is available  
Used for clinical urine catheterization  
Red Latex / Yellow Latex  
Size Length  
Fr8-Fr10 270mm  
Fr12-Fr30 400mm




**GT023-100  
Malecot Catheter**

Made from Natural Latex; Silicone Coated; Length:400mm; Used for Clinical bladder fistula operation urethral catheterization;  
10pcs/box, 200pcs/carton  
Red Latex / Yellow Latex  
FR14-FR32, FR12 and FR34-FR42 are also available for customer requests  
Packing: Individual Blister Bag



**GT022-900  
2-Way All Silicone Foley Catheter**

Made from high quality 100% silicone, medical grade  
With X-ray or without is available  
Fr8-Fr10 (Pediatric)  
Length:310mm Balloon Capacity:3cc/5cc  
Fr12-Fr26 (Standard)  
Length:400mm Balloon Capacity:5cc/10cc  
Fr16-Fr24 (Standard)  
Length:400mm Balloon Capacity:5cc/10cc/30cc  
Packed in individual Blister bag  
10pcs/box, 500pcs/ctn, Dimension:52X43X47.5cm



**GT022-200S  
Super Lubricating Foley Catheter**

Made from Natural Latex; Available With the option of plastic valve or rubber valve; Easy for inserting and hydrophilic coating to reduce the stimulation to mucous membrane of the urethra.  
Type: Single way; 2-way; 3-way; 2-way Female; 2-way Tiemann;  
size: Fr8-24 (Single way)  
Fr8-24 (2-way)  
Fr16-24 (3-way)  
packing: Individual Blister Bag

**GT022-900T  
Disposable All Silicone Foley Catheter with Temperature Sensor**

Made from High quality 100% silicone, medical grade;  
With X-ray or without is available;  
Detect temperature in the bladder,  
temperature range: 25°C-44°C, Accuracy: ±0.1°C  
Length: 400mm; Balloon Capacity: 3ml/5ml/ 10ml/ 15ml/ 30ml  
Fully compatible with most monitors. Temperature extension cable is available.  
Note: Do not pull hard the Temperature measurement cord during using.  
Size: Fr8-24 (Color coded)  
Packing: Individual Blister Bag



**GT022-200T  
Self-Catherization System**

Made from Non-toxic, Non-irritant soft Medical grade PVC; DEHP free is available; Component: Intermittent Catheter; PVC Gloves; Collection bag; Moisture-proof underpad; PVP Swabsticks (3/pkt), sterile.  
Features: Closed system specially designed reducing risk of contamination. The catheter and introducer are coated with a water soluble lubricant for ease of insertion. Catheter Guide allows for simple, protected control of catheter during insertion. 1500ml collection bag eliminates the need for separate container, reducing the chance of spillage.  
Size: Fr6-18 (Color coded)  
Packing: Individual Blister Bag



**GT025-100 T-Drainage Tube**

Made from natural latex, Silicone coated  
Used for the drainage of the operation of Cholecyst  
The connector part is larger in diameter than other parts  
Size: Fr10-Fr26  
Packed in individual Blister Bag



**GT026-100S  
Condom Catheter**

Made from high quality 100% Silicone, medical grade;  
Size: 24mm, 28mm, 31mm, 35mm, 40mm  
With/without adhesive tape  
Packing: Individual Blister Bag



**GT024-100 Penrose Drainage Tube**

Made from natural Latex; Length: 450mm  
Size: 1/4", 5/16", 3/8", 1/2", 5/8", 3/4", 7/8", 1"



**GT026-100 Condom Catherer**

Made from natural latex  
Size: 20mm, 25mm, 30mm, 35mm, 40mm  
Packing: Individual Blister Bag

**GT022-200CP Catheter Plug**

The Spigot is used to provide stop of flow for catheters hygienically during nursing procedures. It is non-invasive which has been utilized to seal the catheter for a short time to allow urine to collect in the bladder.  
• Made from non-toxic PP.  
• Universal size is molded to fit all sizes of catheters or tubes.  
• Graduated design make it to easily accommodate tubing of various internal diameters.  
• Ribbed designed provides strong, secure connections.  
• Ergonomic flanged design assists withdrawal and maximize grip.  
• EO Sterile, Single use only.  
• Total Legth: 55.0±5.0mm Diameter: 12.0±1.0mm

