

DENTAL



**GT081-100
Dental Sponges**

Cotton filled
Cotton Yarn:40's x 40's, 32' x32'
Mesh:19 x 9, 20 x12 etc.
Size:2" x 2", 3" x 3", 4" x 4" etc.
Ply: 8ply etc
Sterile or non-sterile
Packing: 100pcs/ pack, 50pack/ctn
Size: 2" x 2" x 8 ply, 100pcs/pack,
50 pack/ctn, Dim:: 52 x 42 x 45cm;
Size: 3" x 3" x 8 ply, 100pcs/pack,
50 pack/ctn, Dim:: 52 x 50 x 30cm;
Size: 4" x 4" x 8 ply, 100pcs/pack,
50 pack/ctn, Dim:: 52 x 52x 40cm;



**GT080-100
Dental Cotton Rolls**

Description: 100% Cotton,natural,
super Absorbent
A: 20g/bundle, Dia: 7cm
φ 0.8cm X L3.8cm, 56pcs/bundle,
75bundles/ctn
φ 1cm X L3.8cm, 41pcs/bundle,
75bundles/ctn
φ 1.2cm X L3.8cm, 30pcs/bundle,
75bundles/ctn
φ 1.4cm X L3.8cm, 23pcs/bundle,
75bundles/ctn
Packing: 36 X 22 X 20cm
B. 25g/bundle, Dia: 7.2cm
φ 0.8cm X L3.8cm, 60pcs/bundle,
40bundles/box, 12boxes/ctn
φ 1cm X L3.8cm, 50pcs/bundle,
40bundles/box, 12boxes/ctn
φ 1.2cm X L3.8cm, 38pcs/bundle,
40bundles/box, 12boxes/ctn
Packing: 61 X 47 X 43cm



**GT165-100
Salive Ejector**

Material: PVC
Length: 15cm
Packing:100pcs/bag,
10bags/ctn, 31.5X24X16cm



GT094-100 Dental tray

Content Of The Set:
1. Plastic Tray
2. Periodontal Probe With Plastic Handle
3. Choice Of A Plastic-handle Tweezers (stainless Metal Points)
Or A Metal-handle Tweezers Or Forceps.
4. Plastic Stomatoscope, Glass Mirror. Or Pc Mirror, abs Handle.
5. Dental Bib.
6. Additional Supplies Such As A Saliva Ejector, Gel Tray,
Gauze Swabs, Cotton Balls, Dental Brush,
7. Dental Cotton , Mask,etc.
Can Be Assembled To The Individual's Specifications.
8. Packing:1 Set/plastic Bag,200sets/ctn,61x41x39cm,
G.W.:12kg N.W.:10kg
9. Packing:1 Set/blister Bag,480sets/ctn, 61 x 47.5 x 41.5,
G.W.:13kg N.W.:12kg
10. Sterilized With E.o.gas.
11. Expiration Date: 2 Years.



GT094-200A

GT094-200B

GT094-200A / GT094-200B Dental Bib

Item: Disposable Bibs
Material: 2 X 20g Paper Laminated 20g Pe Film
Size: 33 X 46cm, 33 X 66cm Or 40 X 45cm
Color: White /green / Pink/blue/green/purple
Packing: 125pcs/bag,4bags/ctn, 37 X 34.5 X 26.5cm,
G.W.: 4.5kg, N.W.: 4.2kg



Type: Straigh



Type: Curved

**GT038-201
Dental Needle**

Size	18G	19G	20G	21G	22G	23G	24G	25G	27G
Color	Pink or Black	Black	Yellow or Black	Green or Black	Grey or Black	Blue	Purple	Blue or Orange	Grey



**GT038-202
Dental Syringe**

Material: PP
Curved Tip, 12ml



**GT038-200
Dental Needle**

Size: 27g, 30g
Description: Sterilized By Gamma
Packing: 1pc/blister Bag, 100pcs/box, 5000pcs/ctn



**GT038-203
Beauty Syringe**

Material : Medical PC
Luer Lock, 1ml



**GT095-100 Basic Dressing Set 1
(paper-plastics Compound Bag)**

100bags/ctn, 47 X 31 X 41cm
1pc Dressing Tray
5pcs Cotton Wool Balls 05gm
5pcs Gauze Swab 7.5cm X 7.5cm-5pcs
3pcs Dressing Forceps
1pc Surgical Drape
1pc Hand Towel
1pc Limpet Bag
1pc Sterile Field



GT095-200 Basic Dressing Set 2(shrink Bag)

200bags/ctn, 63x29x34cm
Item: 1pc Dressing Tray
8pcs Cotton Wool Balls 0.5gm
5pcs Gauze Swab 7.5cm X 7.5cm-5pcs
3pcs Dressing Forceps
1pc Surgical Drape
1pc Sterile Field



GT095-300 Basic Dressing Set 3(blister Bag)

Item: 1pc Dressing Tray
6pcs Gauze Balls
2pcs Dressing Forceps
1pc Surgical Drape
1pc Hand Towel



GT095-301 Dressing set

120bags/ctn, 41 X 34 X 34cm
1pc Forcep 4#
1pc Surgical blade
3pcs Gauze swab 5cmx5cm



GT095-302 Dressing set

40bags/ctn, 38 X 29 X 28cm
1pc PE+PAPER 33x45cm
2pc Forcep 4#(blue and green)
5pcs 40g/m2(D3cm) 0.8G Non-woven ball



GT095-303 Dressing set

40bags/ctn, 38 X 29 X 28cm
1pc PE+PAPER 33x45cm
2pc Forcep 4#(blue and green)
5pcs 20x20cm(28x24mesh) 1.2G Gauze ball
5pcs Non-woven swab 7.5cmx7.5cm



GT095-304 Dressing set

200bags/ctn, 69 X 42 X 34cm
2pc Forcep 4#(white and green)
3pcs Gauze ball 28x24mesh(15x20cm)



GT095-305 Dressing set

40bags/ctn, 39 X 30 X 30cm
1pc PE+PAPER 32x45cm
1pc Hemostasia Scissor
1pc Forcep 4#
5pcs 20x20cm(28x24mesh) 1.2G Gauze ball
5pcs Non-woven swab 7.5cmx7.5cm



GT095-306 Dressing set

40bags/ctn, 39 X 30 X 30cm
1pc PE+PAPER 33x45cm
1pc Forcep 4#
5pcs 40g/m2(D3cm) 0.8G Non-woven ball



GT124-100 Plastic Forceps

Size: 10.5cm, 12.5cm, 14cm, 15cm, 19cm



GT124-200 Hemostasis Pincers

Size:length 14cm X width 7.5cm
Material:ABS and nylon are available.
ABS for EO sterilization,
nylon material for high temperature steam sterilization.
Color:blue,green,black,red etc. are available.
Packing:individual sterile blister pack