



**GT129-100  
Cervical Scraper (wooden)**

1pc/paper wrap, 100pcs/box  
Non-sterile: 100pcs/box



**GT128-100 Wooden Tongue Depressor  
GT128-200 Plastic Tongue Depressor**

Size Adult: 150mmx18mmx1.6mm  
Sterile: Single Medical Paper Wrapped,  
1pc/paper wrap, 100pcs/box  
Non-sterile: 100pcs/box  
Size Child: 140x14x1.6mm  
Sterile: Single Medical Paper Wrapped,  
1pc/paper wrap, 100pcs/box  
Non-sterile: 100pcs/box



**GT148-203  
Lubricating Jelly Series**


Name	Package	Amount/Unit	Units/Case
Lubricating Jelly 3g	Bag	144/bx	12/cs
Lubricating Jelly 5g	Bag	144/bx	12/cs
Lubricating Jelly 5g	Tube	144/bx	4/cs
Lubricating Jelly 42g	Tube	12/bx	6/cs
Lubricating Jelly 57g	Tube	12/bx	6/cs
Lubricating Jelly 82g	Tube	12/bx	6/cs
Lubricating Jelly 113g	Tube	12/bx	6/cs
Lubricating Jelly 1.0oz	Bag	20/bx	18/cs



**GT113-100 Condom**



Surface: smooth, dotted, ribbed lubricated  
Size: 49mm, 52mm, 55cm  
Be available with spermicide nonoxynol  
Be available with various fruit smell  
If quantity is large enough, the condom is available with various color  
Specification is available with En600, ISO4074, ASTM, WHO, SABS, etc.

# ECG Paper




- Single-lead ECG Paper
- 3-lead ECG Paper
- 6-lead ECG Paper
- 12-lead ECG Paper

**GT148-101  
ECG Paper**




50mmX20m-30m

Single-lead ECG Paper




210mmX300mm-200P

12-lead ECG Paper



6-lead ECG Paper



3-lead ECG Paper

- Single-lead ECG Paper: 50mm X 20m, 200pcs/ctn
- 3-lead ECG Paper: 90mm X 90mm-200P, 100pcs/ctn
- 6-lead ECG Paper: 112mm X 90mm-200P, 100pcs/ctn
- 12-lead ECG Paper: 210mm X 300mm-200P, 20pcs/ctn



## Disposable Electrode

To comply with CE, ISO 13485:2003 and FDA,  
 Comfortable and Convenient, no residual and non-irritable  
 Excellent adhesive from special gel formula  
 Different sizes and material have been studied to meet the different requirements.

### General Monitoring

Foam Recommended for ICU/CCU, Surgery, Emergency			
Item Code	Description	Packaging	Case Dimension
GT148-SF 01	Child, Solid Gel, 36x44mm	50/pk, 40pk/cs	48x36x28cm
GT148-SF 05	Adult, Solid Gel, 42x50mm	50/pk, 40pk/cs	48x36x28cm
GT148-SF 08	Adult, Solid Gel, Ø50mm	50/pk, 40pk/cs	48x36x28cm
GT148-BA 08	With 4mm banana plug, Ø50mm	50/pk, 40pk/cs	43x55x28cm
GT148-SF 07	Neonatal, Solid Gel, Ø25.4mm	50/pk, 40pk/cs	40x36x28cm
Cloth Recommended for ICU/CCU, Emergency, Long term			
Item Code	Description	Packaging	Case Dimension
GT148-SN 08	Cloth, Solid Gel, Ø50mm	50/pk, 40pk/cs	40x36x28cm
Clear Tape Recommended for ICU/CCU, Surgery, Emergency, Long term, Holter			
Item Code	Description	Packaging	Case Dimension
GT148-ST 31	Breathable clear tape, Solid Gel, Ø50mm	50/pk, 40pk/cs	40x36x28cm

### Radio Translucent

Foam Recommended for Surgery, Emergency, Cath Lab, Radiology			
Item Code	Description	Packaging	Case Dimension
GT148-SC 08	Adult, Solid Gel, Ø50mm	50/pk, 40pk/cs	40x36x28cm

### Stress / Diaphoretic

Recommended for ICU/CCU, Surgery, Stress, Holter, Emergency			
Item Code	Description	Packaging	Case Dimension
GT148-SF 10	Adult, Solid Gel, 42x37mm	50/pk, 40pk/cs	40x36x28cm

### Wet Gel


Recommended usage ICU/CCU, Surgery, Stress, Holter, EMS, Emergency			
Item Code	Description	Packaging	Case Dimension
GT148-WF 32	Adult, Wet Gel, 52x47mm	50/pk, 40pk/cs	40x36x28cm

### Repositionable

Recommended usage ICU/CCU, Surgery, Emergency, Cath Lab, Radiology			
Item Code	Description	Packaging	Case Dimension
GT148-RS 01	Adult, Solid Gel, 32x40mm	3/pk, 10pk/bx, 12bx/cs	48x22x22cm

## Disposable Electrosurgical Pencil

- 170mm-ABS Pencil appearance
- 7cm Stainless steel tip
- 2.36mm-diameter of the tip
- 3m-PVC medical cable
- Maximal out put current:1A
- Frequency confine:0-1.0MHZ
- Ways of aseptic technique: Sterile EO
- Current resistance:200



### GT148-200 Medical Conductive Paste

Property Unit Value  
 Sound Speed(35°C)M/s1520-1620  
 Electric Resistance(35°C)-13-17  
 Unit Area Acoustic Impedance (35°C)  
 Pa-s/m1.5-1.7x106 Sound Attenuation  
 (35°C)Db/(cm.mhz)≤0.05  
 Viscosity (25°C)Pa.s≥15  
 Ph-5.5-8.0  
 Size: 60ml, 100ml, 150ml, 250ml



### GT148-201/202 Ultrasonic Gel

Property Unit Value  
 Sound Speed(35°C)m/s1520-1620  
 Unit area acoustic impedance (35°C)  
 Pa-s/m1.5-1.7x106 Sound attenuation  
 (35°C)dB/(cm.Mhz)≤0.05  
 Viscosity (25°C)Pa.s≥15  
 PH-5.5-5.8  
 Size: 250ml, 1000ml, 2500ml, 5000ml