

Homecare



GT103-100 Hot Water Bag

Size: 750cc/ 1000cc/ 1750cc/ 2000cc
/ 2500cc/ 3000cc
Color: Blue, Green, Red
Material: PVC
Double Sided Trill, Single Sided Trill
Packing: 1pc/ 1 Pe Bag



GT103-100 Hot water bag

Size: 750cc/ 1000cc/ 1750cc/ 2000cc/ 2500cc/ 3000cc
Color: Blue, Green, Red
Material: Latex ≥ 30%,
Bs Standard: latex ≥ 50%,
Double Sided Trill, Single Sided Trill
Packing: 1pc/ 1 Pe Bag



Gt103-200 Cover For Hot Water Bag

Material: Flannelette/ Knit/ Furry Etc.
Size: 750cc. 1000cc. 1500cc. 1750cc. 2000cc. 2500cc
Packing: Individual Opp Bag



GT104-100 Ice Bag

Size: 6"/ 9"/ 11"
Packing: PVC Box/ PE Bag/ Paper Box
Can be customized design.



GT115-001 Instant Ice Pack

280G, PE+PA, 23x14CM
160G, PE+PA, 19.5x12CM
100G, PE+PA, 15x11CM

Usage:

Find the liquid bag in the ice bag and make it broken and 3 seconds later it will be frozen. Then shake the bag up and down to make the content fully mixed, and the temperature will be depressed below zero by the ice bag in 2 minutes.

Purpose:

To control and alleviate the gore and pain caused by strain, bruise, pull and burn, also can rapidly remove the pain and discomfort induced by fever, headache, toothache, mosquito bite.

Freeze and keep fresh vegetable, fruit, beer, beverage and food and so on when no cold source available.

Precaution:

The ice bag can only be used for external medication. If eye or skin should contact the content, it must be washed by clean water thoroughly. If one eats the content by mistake, he must drink enough water, try his best to vomit it and turn to doctor when necessary.

In order to avoid a too low temperature, it should be better to wrap the ice bag with a towel or cotton cloth to begin with. Those who have something wrong with their circulation system should refer to doctor in advance.

Characteristic:

This product can refrigerate instantly without cold source and be easily used because of portability.

GT115-100**Instant Hot & Cold Pack**

Material: PVC+Sodium Acetate+Sheet iron
 Size: 9x9cm 61g; 10x11.5cm 115g; 11.5x11cm 85g;
 14x13.5cm 185g; 9.5x10.5cm 90g
 Packing: Individual OPP bag

**GT115-100****GT115-101****GT115-101****Instant Hot & Cold Pack**

Material: PVC+Liquid Gel
 Size: 9.5x8cm 40g; 13x13cm 100g; 10x11cm 60g;
 17x13cm 200g; 10x10cm 65g
 Packing: Individual OPP bag

GT115-105**Instant Hot & Cold Pack**

Belt
 Material: PVC+Round Gel ball+Nylon
 Size: 27x31; 13.5x33"

**GT105-100****Cold Eye Mask**

Material: Outer bag: PVC
 Infilling: SAP Ice ball+glycerinum+Distilled water
 30G, 20.5x5.2CM
 120G, 77x9CM
 150G, 68x10CM
 55g, 24x7.5cm
 56g, 65x8cm

**GT115-102/104****Instant Hot & Cold Pack**

Material: Nylon Taffeta+PVC
 250G, 23x13cm
 200G, 20.5x15.5cm

**GT105-100 Air Cushion**

Item: Air Cushion
 Material: Calcium carbonate + 30% rubber
 Description:

| | | | |
|----------|---------------|-----------|------------|
| φ 35cm | 1pc/ 1 PE bag | 20pcs/ctn | 37X37X13cm |
| φ 40cm | 1pc/ 1 PE bag | 20pcs/ctn | 42X42X15cm |
| φ 42.5cm | 1pc/ 1 PE bag | 20pcs/ctn | 45X45X16cm |
| φ 45cm | 1pc/ 1 PE bag | 20pcs/ctn | 47X47X17cm |

**GT106-100 Ear Syringe**

Item: Ear Syringe
 Material: Rubber/PVC
 Description:

| | | | |
|------|---------------|------------|------------|
| 30ml | 1pc/ 1 PE bag | 600pcs/ctn | 37X34X40cm |
| 60ml | 1pc/ 1 PE bag | 600pcs/ctn | 64X34X40cm |
| 90ml | 1pc/ 1 PE bag | 360pcs/ctn | 55X34X40cm |

**GT107-100 Breast Reliever**

Item: Breast Reliever
 Material: PVC + 25% rubber
 Packing: 1pc/1 bag,
 144pcs/ctn, 63X47X48 cm

**GT108-100 Vaginal Syringe**

Material: PVC
 Size: 90g/150g/220g

**GT109-100 Pear Shape Syringe**

Material: PVC
 Size: 300g



Plastic Pill Box



GT164-100
Lover Shape 4 Case Pill Case
Material: PP
Size: 8.3x8.3x2.1cm




GT164-101
4 Case Pill Case
Material: PP
Size: 12x3x3.2CM



GT164-102
8 Case Pill Box
Material: PP
Size: 10x7x2cm



GT164-103
14 Case Pill Box
Material: PP
size: 18.3x9.3x3.3cm



GT164-104
7 Case Spring Pill Orgnaizer
Material: ABS+PP
Size: 19.7x4.7x2.3cm



GT164-105
14 Case Detachable
Material: PP
size: 17x7.2x2.5cm



GT164-106
31 Case Pill Box
Material: HIPS+PP
Size: 28.5x10.5x2.8cm



GT164-107
28 Case Pill Box
Material: ABS+PP+PS
Size: 21x12.5x2.3cm



GT164-108
6 Case Pill Box
Material: PP
Size: 9x8x3cm



GT164-109
The Detachable 7 Case Pill Box
Material: PP
size: 17.5x3.5x2.5cm



GT164-110
14 Case Pill Box
Material: PP
Size: 13.5x7x2.5cm



GT164-111
4 Case Pill Box
Material: PP
Size: 9.5x3.5x2.5cm



GT164-200
Pill Crusher
Material: POM+ PS
Feature Crush And Store Pills



GT164-201
Pill Crusher
Material ABS+ PS
Feature Crush And Store Pills
Colour Blue; Red; White
Function Easy To Use For Pill
Storage And Crushing



GT164-202
**Multifunction Medical
Pill Crusher**
Material: ABS+Stainless Steel Blade



GT164-112
Keyboard Pill Box
Material: PP
size: 23x13x3.2cm



GT164-113
28Case Pill Box
Material: PP
Size: 20x11x3.8cm



GT164-114
Oval 28Case Pill Oegnaizer
Material: ABS+PP+PS
size: 9x9x10cm



GT164-115
28Case Pill Oegnaizer
Material: PP +ABS
size: 9.6x5.9x12cm



GT164-116
Drawer 28Case Pill Box
Material: ABS+PS
size: 10.5x12.6x4.3cm



GT164-203
Pill Cutter
Material: PS
Size: 9.5x4.5x3cm



GT164-204
**Multi-function Group
of Storage Bottles**
Material: Tritan+HIPE
Size: 6x6x20cm