

CONTENT

Company Profile 02-05	Hospital Bed 06-22	Transport Stretcher 23-26
Operating Room Equipment 27-39	Gynecology Equipment 40-53	Medical Trolley 54-56
	Other Products 57/60	
	After Sale Service 61	



To be the best hospital furniture supplier

WE ARE DIFFERENT

CHANGE and REBIRTH

TO BE THE BEST HOSPITAL FURNITURE SUPPLIER

NO SURPRISE NO SALE

KEEP MOVING ALL IN

www.spectrum-medlab.com



AG-BR002C

7功能ICU病床
7-Function ICU Bed

Specification:

1. Electrostatic powder coating for the whole surface, 10-part cold rolled steel plate punching bed board, Cold-rolled steel bed frame.
2. Linak Motor from Denmark, three column straight lift system
3. Alarm system
4. 4-piece ABS big handrails, ABS detachable headboard and footboard, embedded handrail operator.
5. Electronic weighing system accurate to 0.01kg, with memory function, operate easily
6. CPR board quick release
7. Continue working for 8 hours without battery support
8. Foot-operated brake control system, four 125mm double silent castors
9. ABS shell one time formed dust cover, angle indicator

1. One-button electric CPR
2. One-button reset
3. One-button chair (similar to TV-position)
4. Night light
5. With weighing scale





AG-BY003C

5功能电动病床
5-Function Electric ICU Bed

1. Backrest max upward angle: 75°
2. Footrest max upward angle: 45°
3. Height adjust: 480-760mm
4. Trendelenburg: 12°
5. Anti-trendelenburg: 12°

尺寸 Size :
2200 × 1100 × 480-780mm



CPR功能
CPR function, accessible from both side of the bed




双面静音脚轮
Double-side silent wheels with locking



可拆卸ABS床板
4-Part of detachable ABS bedboard



进口丹麦电机
Three years warranty for Linak motor



不锈钢紧固件
All screws are made of stainless steel




酒红色 Wine red



蓝色 Blue



木纹色 Wooden



AG-BR001

8功能电动病床
8-Function ICU Bed

尺寸 Size :
2000-2200 × 1020 × 525-965mm



- | | |
|--------------------------------|----------------|
| 8-功能 8-Functions: | |
| 高度调节 HI-LOW | 525-965mm |
| 背板角度调节 Back section adjustable | 0-75° |
| 腿板角度调节 Knee section adjustable | 0-45° |
| 前倾 Trendelenburg | 0-15° |
| 后倾 Anti-Trendelenburg | 0-15° |
| 床面延展 Footboard extend function | 2000-2200 |
| 座椅功能 Chair function | 0-30° |
| CPR 功能 CPR function | |
| 选配 Optional Accessory | Weighing scale |



AG-BY009

5功能电动病床 (称重型)
5-Function ICU Bed (Weighing Type)

尺寸 Size :
2250 × 1050 × 575-815mm



AG-BY009重症监护病床是专门为重症患者设计的具有称重系统精确到0.01公斤功能的多功能电动病床,配有10片式床板,舒适性更佳。嵌入式操作板设计使病人和护士操作更加方便。爱琴海专利设计插头,能够适应各种插座并且适用于多种电源电压。双立柱结构垂直升降无噪音,电机使用寿命更长久。

AG-BY009 ICU bed is designed for intensive care patients with weighing function accurate to 0.01kg. Equipped with 10pcs bedboards, more comfortable for patients to lay on. Embedded operator design more convenient for patients and nurses to operate. AEGEAN patent plug design, can replace all the plug and apply to any voltage. Vertical lift noiseless, extend the life of motor.

- | | |
|--------------------------------|-------------------------|
| 5-function: Weighing function | |
| Backrest max upward angle: 75° | Trendelenburg: 12° |
| Footrest max upward angle: 40° | Anti-trendelenburg: 12° |
| Height adjustment: 575-815mm | |



AG-BY007

5功能电动病床
5-function Electric Bed

- | | |
|-------------|--------------------------------|
| 尺寸 Size : | 2150 × 1090 × 450-720mm |
| 5-function: | Backrest max upward angle: 75° |
| | Footrest max upward angle: 40° |
| | Height adjustment: 450-720mm |
| | Trendelenburg: 12° |
| | Anti-trendelenburg: 12° |



AG-BM005

5功能电动病床
5-function Electric Bed

- | | |
|-------------|--------------------------------|
| 尺寸 Size : | 2000 × 970 × 500-750mm |
| 5-function: | Backrest max upward angle: 75° |
| | Footrest max upward angle: 40° |
| | Height adjustment: 500-750mm |
| | Trendelenburg: 12° |
| | Anti-trendelenburg: 12° |



AG-BY004

5功能电动病床
5-Function ICU Bed

- | | |
|-------------|--------------------------------|
| 尺寸 Size : | 2200 × 1060 × 440-760mm |
| 5-function: | Backrest max upward angle: 75° |
| | Footrest max upward angle: 40° |
| | Height adjustment: 440-760mm |
| | Trendelenburg: 12° |
| | Anti-trendelenburg: 12° |



AG-MC001

5功能电动护理床 5-Function Electric Nursing Bed

背部调节, 腿部调节, 前后倾斜, 升降调节
可拆卸折叠, 方便贮存

5寸脚轮

尺寸: 2150 × 1030 × 460-860mm

承重170kg

进口林耐克电机



Back section adjustable

Foot section adjustable

High-low 460-860mm

Trendelenburg, anti-trendelenburg

1. Depatchable, foldable. save space when shipping and stocking

2. 125mm dia castors

3. Size: 2150*1030*460-860mm

4. Loading:170kg

5. Four Linak motors from Denmark.



AG-W001

超低位五功能电动护理床 Ultra-low Electric Homecare Bed

可选配件

1、吊臂

2、配套床垫可选: 5mm棕垫+充气防褥床垫

3、护士操作器(带电动CPR)

4、配套床头柜与床边桌

5、可选4寸中控刹车脚轮或者可调节脚垫

6、背板可选透X光电木板背板

7、分段式护栏可选往中间倒、往两边倒、或者六栏式一体护栏

Optional accessories:

1. Lifting pole

2. Optional mattress:5mm coir mattress+inflatable mattress

3. Nurse operator(with electric CPR function)

4. Matching cabinet and over bed table

5. 4-inch central control braking castors or adjustable foot pad

6. Backboard plate with X-ray permeable

7. Segment guardrail in the middle or integrated six bars guardrail

尺寸 Size:

2200 × 1100 × 250-650mm



AG-WS001

三功能手动护理床 Three-Function Manual Homecare Bed

尺寸 Size:

2150 × 1100 × 350-600mm



颜色可选
Colors optional



AG-BMY002

三功能液压床
3-Function Hydraulic Bed



尺寸 Size:
2160 × 970 × (450 - 700)mm
Backrest controlled by gas spring : 0-75°
Footrest controlled by gas spring : 0-40°
Height adjustable controlled by hydraulic : 450-700mm

Technical parameters:

1. Frame : cold-rolled steel structure, electro-coating and powder-coating
2. Four - piece carbon steel stamping bed, cold-rolled steel frame
3. American power hydraulic cylinder
4. 4 pieces of ABS side handrails, ABS detachable head/foot boards
5. Dynamic loading: 250KG, Static loading: > 300KG



AG-BY104

三功能手电一体病床
3-Function Manual
and Electric Bed

尺寸 Size:
2200 × 1060 × (440-760)mm
Backrest max upward angle: 75°
Footrest max upward angle: 40°
Height adjustment: 440-760mm





AG-BM005

5功能电动病床
5-Function Electric Hospital Bed

尺寸 Size:
2000x970x(500-750) mm



AG-BMS002B

3功能手摇床
3-function manual bed(ABS bedboard)

背部抬起0-75°
腿部抬起0-40°
整体上升下降460-760mm

材料说明:

床体表层静电粉末喷涂
四片式ABS床板
六档铝合金护栏, ABS可拆卸床头尾板
304不锈钢可调节输液杆
4个直径125mm双面静音脚轮
所有螺丝均为不锈钢材质
脚踏式中控刹车系统, 引流挂钩
动态承重: 250Kg
静态承重: >300Kg

Backrest adjustment:0-75°
Footrest adjustment:0-40°
HI-LOW:460-760mm

Technical parameters:

- 1.Treated by electro-coating and powder-coating.
- 2.4-part ABS bed boards.
- 3.6-crank al-alloy handrails, ABS head/foot board.
- 4.Height adjustable 304 stainless steel IV pole.
- 5.4 double-side silent wheels(125mm)
- 6.All the screws made of stainless steel material.

Dynamic load: 250Kg
Static load:> 300Kg



AG-BMS001

5功能手动病床
5-Function Manual Bed

尺寸 Size:
2160 x 960 x (500-750)mm





MANUAL BED

AG-BYS106

2功能手动病床
2-Function Manual Bed

尺寸 Size:
2120 × 1050 × 540mm



AG-BMS101A

2功能手动病床
2-Function Manual Bed

尺寸(Size):
2130 × 970 × 540mm



ISO 13485:2003 CE

AG-BYS203

1功能手动病床
Single Function Manual Bed (Type I)

尺寸 Size:
2130 × 960 × 540mm

功能 Functions:
背板角度调节
Back section adjustable: 0-75°



AG-BYS204

1功能手动病床
Single Function Manual Bed (Type II)

尺寸 Size:
2130 × 960 × 540mm

功能 Functions:
背板角度调节
Back section adjustable: 0-75°



AG-BMS303

平床
Flat Bed

尺寸 Size:
2130 × 960 × 540mm





AG-CB001

1功能手动病床
1-Crank Manual Child bed

技术参数

1. 尺寸: 1800 × 820 × 500mm
2. 耐用框架: 静电喷涂和粉末涂层处理的冷轧钢

功能

背板可调: 0-75°

Technical Parameter

1. Overall size: 1800 × 820 × 500mm
2. Durable frame: made of cold-rolled steel plate, treated by electro-coating and powder-coating

1-function:

Backrest max upward angle: 0-75°



AG-CB004

儿童床
Flat Child Bed

- 1) Overall size (L*W*H): 1210*690*600mm
- 2) Durable frame: made of cold-rolled steel plate, treated by electro-coating and powder-coating
- 3) Detachable Steel headboard/footboard with plastic board
- 4) 2-part Steel bedboards
- 5) Full length steel handrails
- 6) Fully tested sturdy construction capable of taking a maximum user weight of up to 200kg
- 7) With one ss height adjustable I.V pole, with Drainage hook
- 8) Four silent wheels with cross brakes



AG-CB003

儿童床
Flat Child Bed

技术参数

1. 尺寸: 1810 × 900 × 590mm
2. 耐用框架: 静电喷涂和粉末涂层处理的冷轧钢

Technical Parameter

1. Overall size: 1810 × 900 × 590mm
2. Durable frame: made of cold-rolled steel plate, treated by electro-coating and powder-coating



AG-CB014

2功能婴儿床
2- Function Child Bed

Technical Parameter

1. Size: 800 × 700 × 1600mm
2. 304 Stainless steel bed frame
3. Footboard with plastic patient name card
4. 4-part cold-rolled steel bedboards
5. Detachable 304 stainless steel handrails
6. Max weight 250Kg
7. Castors with cross braking system
8. Stainless steel retractable IV pole
9. Drainage hook



AG-CB009

医用婴儿车
Baby Crib(Type II)

尺寸 Size:
810 × 520 × 940mm

Tilting angle of head from 0 to 12 °

Material:

Steel bed frames with basket
Steel side rails around
ABS tray with mattress
4 castors

Maximum loading capacity: 60kgs



AG-CB011

医用婴儿车
Baby crib

尺寸 Size:
850 × 500 × 780-980mm

Technical Parameter


1. Steel frame
2. ABS baby basin
3. ABS tray and base
4. 4 castors with cross brakes
5. Tilting angle of head from 0 to 12 °
6. Height of bed adjustable from 810mm to 930mm
7. Maximum loading capacity: 60kgs

1. All Aegean hospital beds use Linak Motor.
电机可选 Motor Optional: Linak ,LK ,JC,T-motion
- 2.Package: wooden or carton
- 3.OEM medical is available
- 4.Electric Bed 20ft:16pcs 40HQ:55pcs

5. 颜色可选: Color Optional:




防褥疮气垫
AG-M001
Anti-decubitus Air Mattress(Bubble)



防褥疮气垫
AG-M002
Anti-decubitus Air Mattress(Stripe)



医用床垫
AG-M004
4-Folding Foam Mattress



医用床垫
AG-M011
Medical Mattress



不锈钢静脉输液架
AG-IVP004 Stainless Steel IV pole




输液架
AG-SS009 IV Stand



输液架
AG-P027 IV Pole



吊臂
AG-P020 Lifting Pole



10cm



ISO 13485:2003 CE

AG-HS001

液压担架车

Hydraulic transport stretcher

尺寸 Size:
2100 × 850 × 600-950mm

材质:
液压泵: 美国进口液压泵
钢制框架, 铝合金护栏
背板可X-ray
底部罩壳为进口PP材质
配IV Pole

Material:
Pump: American Pump
Steel frame, aluminum alloy handrails
X-ray available back board
Imported PP material bottom platform
With I.V. pole



AG-HS010

电动担架车
Electric Transport Stretcher

尺寸 Size:
1930 × 700 × (630-930)mm



AG-HS021

8轮手术交换车
Exchange Stretcher Trolley For Operation Room

功能:

最大程度代替了手术室对接车, 速度更快, 一车完成所有繁琐的环节, 8轮设计, 避免灰尘及细菌接触手术室地面

Functions:

- Back section adjustable: 0-90°
- Height adjustable
- 8 wheels exchange for OT Room

尺寸 Size:

1900 × 640 × 730-940mm



AG-HS003

手术交换车
Conneting Stretcher Trolley
For Operation Room



技术参数:

配PP护栏, 输液架

Technical Parameter:
PP handrails, IV pole

功能 Functions:

- HI-LOW: 640-970MM
- Back section adjustable: 0-75°

尺寸 Size:

3650 × 640 × 640-970mm





AG-HS002
手动担架车
Manual Transfer Stretcher

尺寸(Size):
1930 × 640 × (540-840)mm



AG-HS008
液压担架车
Hydraulic Stretcher

尺寸(Size):
1930 × 700 × 630/900mm



AG-HS013
手动担架车
Manual Transport Stretcher

尺寸(Size):
1900 × 630 × 800mm



AG-HS014
手动担架车
Manual Transport Stretcher

尺寸(Size):
1900 × 630 × 600-800mm



AG-OT008
DR专业手术台
DR Professional Operating Table

主要功能:
电动液压控制, 四点式自我补偿式刹车系统, 平移530mm, 升降450mm。床面可形成V形体位100度, 反V形体位220度。

Main Usage:
Electric-hydraulic control four point self compensation.

台面尺寸 Dimension of table top:
长 Overall Length 2100mm
宽 Width 500mm

部件尺寸 Dimension of section:
头板尺寸 Head section 250 × 500mm
背板尺寸 Back section 465 × 500mm
腰桥 Waist section 75 × 500mm
坐板尺寸 Seat section 630 × 500mm
腿板 Leg section 600 × 240mm

主要技术参数
台面最低 Min height 700mm
台面最高 Max height 1150mm
升降范围 Adjustment 450mm
前倾 Reversed trendeleburg $\geq 30^\circ$
后倾 Trendeleburg $\geq 25^\circ$
左右倾 Lateral tilt $\geq 25^\circ$
前后平移 Move forward and backward ≥ 530 mm
头板上折 Head section raised from the horizontal $\geq 45^\circ$
头板下折 Head section lowered from the horizontal $\geq 90^\circ$
背板上折 Back section raised from the horizontal $\geq 75^\circ$
背板下折 Back section lowered from the horizontal $\geq 30^\circ$
腿板下折 Leg Section lowered from the horizontal $\geq 90^\circ$
腿板外展 Leg section opened $\geq 90^\circ$
腰桥升降 Kidney bridge elevated ≥ 80 mm
升降速度 Lifting speed 4 ~ 16mm/s



搁腿架
Leg support



搁手板
Arm support



缚身带
Strap



麻醉帘
Anaesthesia screen frame



腰架
Waist fixer



肩架
Shoulder fixer

AG-OT007

高级电动液压手术台
综合豪华手术台 眼科手术
Electric-Hydraulic Operating Table

主要用途:
电动手术台供临床机构手术时使用。

Main Usage:
It can be operated for ophthalmology, gynecology orthopedics clinical all kinds of surgeries.

结构性能:
手术台由机械及电动液压控制系统组成。该产品通过人工按钮操纵，能够满足临床机构手术时所需体位的要求。

Structure :
Operating table consists of two parts: mechanical and electric-hydraulic control system. Depending on manual operation, the product can meet the requirements of posture in clinical operation.

台面尺寸	Dimension of table top	
全长	Overall length	2050mm
宽	Width	550mm

台面分段尺寸	Dimension of section	
头板	Head section	250 × 500mm
背板	Back section	465 × 500mm
腰板	Waist section	75 × 500mm
座板	Seat section	550 × 500mm
腿板	Leg section	600 × 240mm

台面高度	Table top from floor	
台面离地最低	Min height	580mm
台面离地最高	Max height	930mm
行程	Adjustment	350mm



台面调节范围 Range of movements

前倾	Reversed trendelenburg	≥20°
后倾	Trendelenburg	≥15°
左右倾	Lateral tilt	≥15°

台面分段调节最大范围 Adjustment of sections

头板上折	Head section raised from the horizontal	≥45°
头板下折	Head section lowered from the horizontal	≥90°
背板上折	Back section raised from the horizontal	≥75°
背板下折	Back section lowered from the horizontal	≥30°
腿板下折	Leg section lowered from the horizontal	≥90°
腿板外展	Leg section opened	≥90°
腰板升高(选配)	Kidney bridge elevated(Optional)	≥80mm
床台水平旋转	Bed pan rotation	0 ~ 360°
水平复零(选配)	Restoring horizontal position(optional)	±2°
升降速度	Lifting speed	4 ~ 16mm/s
安全工作载荷	Safe working load	250kg



AG-OT009

高级电动手术台
泌尿科手术台 综合手术台
Electric Operating Table (C-arm)



参数 Specifications:

台面的尺寸 Dimension of table top:

长 Overall length:	2100mm
宽 Width:	500mm
头板 Head section:	250 × 500mm
背板 Back section	545 × 500mm
坐板 Seat section	450 × 500mm
腿板 Leg section	600 × 240mm
腰桥 Kidney bridge	75 × 500mm
脚板 Foot section	180 × 273mm
后倾 Reversed trendelenburg	≥ 20°
前倾 Trendelenburg	≥ 15°
横向倾斜 Lateral tilt	≥ 15°
向后和向前运动 Backward and forward movements	≥ 400mm

调整部分

头板角度调节 Head section raised from the horizontal	≥ 45° — ≥ 90°
背板角度调节 Back section raised from the horizontal	≥ 45°
腰桥 Kidney bridge elevated	≥ 80mm
腿板下降 Leg section lowered from the horizontal	≥ 90°
腿部外展 Leg section opened	≥ 90°

AG-OT012

电动手术台
Electric Operating Table



参数 Specifications:

台面的尺寸 Dimension of table top:

长 Overall length:	2100mm
宽 Width:	480mm

部件尺寸 Dimension of section

头板 Head section:	280 × 325mm
背板 Back section	500 × 480mm
座板 Seat section	530 × 480mm
腿板 Leg section	510 × 236mm
脚板 Foot section	180 × 273mm
腰桥 Kidney bridge	100 × 480mm

从台面至地板距离 Table top from floor

最小高度 Min height	750mm
最大高度 Max height	1000mm
调整 Adjustment	250mm

台面活动范围 Range of movements

台面前倾角度 Reversed trendelenburg	≥ 20°
台面后倾角度 Trendelenburg	≥ 15°
台面横向倾角度 Lateral tilt	≥ 15°

头板转折角度 Head section raised from the horizontal	≥ 45°
头板后折角度 Head section lowered from the horizontal	≥ 90°
背板转折角度 Back section raised from the horizontal	≥ 75°
背板后折角度 Back section lowered from the horizontal	≥ 30°
腿板下折角度 Leg section lowered from the horizontal	≥ 90°
腿板打开角度 Leg section opened	≥ 90°
腰桥 Kidney bridge elevated	≥ 100mm

AG-OT015

手动液压手术台
Manual Hydraulic Operating Table (C-arm)



参数 Specifications:

台面尺寸 Dimension of table top:

长 Overall length:	2100mm
宽 Width:	500mm
头板 Head section:	250 × 500mm
背板 Back section	465 × 500mm
腰桥 Kidney bridge	75 × 500mm
座板 Seat section	450 × 500mm
腿板 Leg section	600 × 240mm
脚蹬 Foot rest	273 × 180mm
后倾 Reversed trendelenburg	≥ 15°
前倾 Trendelenburg	≥ 20°
横向倾斜 Lateral tilt	≥ 15°
头板 Head section:	≥ 45° — ≥ 90°
背板 Back section	≥ 45°
腿板垂直 Leg section—lowered from the horizontal	≥ 90°
腿板外摆 Leg section—opened	≥ 90°
前移后移 Backward and forward movements	≥ 400mm
腰桥 Kidney bridge elevated	≥ 100mm
承重量 Loading capacity	250kg

AG-OT025

手动液压手术台
Manual Hydraulic Operating Table



参数 Specifications:


长 Overall length:	2050 ± 50mm
宽 Width:	480 ± 20mm
升降高度 Max & min height of the platform	(750 ~ 950) ± 50mm
背板向上折叠 Back board fold upward	≥ 80°
背板向下折叠 Back board fold downward	≥ 10°
前倾 Reversed trendelenburg	≥ 25°
后倾 Trendelenburg	≥ 25°
左右倾 Leftward and rightward	≥ 20°
头板向上折叠 Head board fold upward	≥ 45°
头板折叠向下 Head board fold downward	≥ 90°
腿板折叠向下 Leg board fold downward	≥ 90°
腿板可分开 Leg of the board split	≥ 180° (可拆卸 dischargeable)
腰板可移动 Waist board lifting	≥ 120mm



AG-45

多臂吊塔
Double Arm Electric Surgical Pendant

- Dry pendant(1pc):
- 1.Horizontal rotation angle $\leq 340^\circ$
 - 2.With three layers instrument platform(height can be adjusted)
 - 3.Standard accessory for the gas connecting: Oxygen(1pc), VAC(1pc), Air(1pc)
 - 4.Electrical outlet: 6 pcs(220V 10A)
 - 5.With gas braking system
 - 6.With equipotential ground terminal(2pcs), one drawer
 - 7.With one adjustable stainless steel IV pole and adjusted stainless steel injection pump shelf
- Wet pendant(1pc):
- 1.Horizontal rotation angle $\leq 340^\circ$
 - 2.With two layers instrument platform(height can be adjusted)
 - 3.Standard accessory for the gas connecting: Oxygen(1pc), VAC(1pc), Air(1pc)
 - 4.Electrical outlet: 6 pcs(220V 10A)
 - 5.With equipotential ground terminal(2pcs)
 - 6.With one adjustable stainless steel IV pole and adjustable stainless steel injection pump shelf



AG-40D

多臂吊塔
Double Arm Electric Surgery Pendant

- 1.Workbench vertical range is 300 mm;Electric power ≤ 1 kw
- 2.Horizontal rotation angle $\leq 340^\circ$;
- 3.Gas terminal standard: Oxygen(1pc),Vacumn(1pc), Air(1pc)
- 4.electric socket: 8pcs,grounding terminal:2pcs
- 5.Adjustable stainless steel iv pole and syringe pump support 1set



AG-360S

单臂电动吊塔
Single Arm Electric Hospital Pendant



AG-40H-1S

单臂电动腹腔镜塔
Single Arm Electric Abdominal Cavity Pendant

- Specification:
- 1 Working platform vertical lift range: 600mm
 - 2 Horizontal rotation angle $\leq 340^\circ$,with gas braking system
 - 3 Loading capacity:100~200KGS
 - 4 Standard accessory for the gas connecting:Oxygen(2pcs), VAC(2pcs), Air(1pc), CO2(1pc)
 - 5 Electrical outlet:8 pcs(220V 10A)
 - 6 With equipotential ground terminal(2pcs), Internet connector RJ45(1pc)
 - 8 Surface treatment using high-pressure electrostatic powder coating
 - 9 Adjustable stainless steel iv pole and syringe pump support 1set

AG-58

旋转吊塔
Surgical Rotary Pendant

- 1 Work bench vertical lift 300MM,motor power ≤1kw,Cross arm lifting,speed 22mm/s
- 2 Radius of the cross arm ≤1000mm, Horizontal rotation angle≤340°
- 3 Gas interface standard configuration: 2 oxygen gas, 2 VAC
- 4 Electrical outlet:8 pcs
- 5 equipotential ground terminal 2pcs,net work port RJ45
- 6 Suspension anesthesia machine
- 7 Surface treatment using high voltage electrostatic powder coating
- 8 Adjustable stainless steel iv pole and syringe pump support 1set



AG-350

单臂电动麻醉塔
Single Arm Electric Anesthesia Pendant

- 1 One pc rotating horizontal arm(length more than 800mm); Horizontal rotation angle≤340°
- 2 Standard accessory for the gas connecting: Oxygen 2,Vacumn2
- 3 Electrical outlet:6 pcs
- 4 With gas braking system, with one drawer, with equipotential ground terminal(2pcs)
- 5 Adjustable stainless steel iv pole and syringe pump support 1set



AG-20V-1

多功能医疗柱
Multi-Function Medical Column

- 1 loading capacity:200KG
- 2 layers: 3pcs(height adjustable); round corner anti-collision design
- 3 Gas outlet standard: 2 Oxygen 2 vacuum 2 Air
- 4 Electric sockets:12
- 5 ground terminal: 4 / internet port:1
- 6 Adjustable stainless steel iv pole and syringe pump support 1set
- 7 1 drawer and 2 basket and 2 lighting lamp



AG-HBD002

医用设备带
Bed Head Unit

- 尺寸 Size:
2000mm × 270mm × 70mm
Head Unit (luxurious) in gas-electric separation design



国标 JIS



德标 DIN



氧气吸入器 Oxygen inspirator



英标 BS



美标 ANSI



吸引缓冲瓶 Vacuum buffer bottle

AG-LT014B

LED无影灯 (双头)
LED shadowless operation (lamp-double)

Illuminance:LED3	≥ 120,000Lux
Illuminance:LED5	≥ 160,000Lux
Color temperature	4500 ± 200K
Colour reduction index(Ra)	96
Illuminance depth	1020mm/1150mm
Rated power of bulb	3.2V/1W
Total irradiance:LED3	380W/m ²
Total irradiance:LED5	510W/m ²
Size of light field	200 ~ 350mm
Service life of illuminant	50,000h
Power supply voltage	AC110-240v,50/60Hz
Lowest height of installation	2900mm
Total power consumption	140W
Total LED bulb quantity	294 pcs



AG-LT019 (Double 75/55)

LED手术灯(双头75/55)
Operation Lamp

Reflector diameter	750mm	550mm
Center illuminance	40~160klx	40~140klx
Spot diameter(d10)	230mm	150mm
Spot diameter(d50)	160mm	120mm
Light beam depth	700mm	700mm
Color temperature	4200 ± 300k	4200 ± 300k
Color rendering index	Ra>95	Ra>95
Maximum illuminance	Ec<1000W/m2	Ec<1000W/m2

Optional Accessory: Camera, Arm with screen



AG-LT019A(Single 55)

LED手术灯(单头55)
Operation Lamp

Reflector diameter	550mm
Center illuminance	40~140klx
Spot diameter(d10)	150mm
Spot diameter(d50)	120mm
Light beam depth	700mm
Color temperature	4200 ± 300k
Color rendering index	Ra>95
Maximum illuminance	Ec<1000W/m2



AG-LT010-1

立式冷光五孔手术无影灯
Operating Lamp

采用从不同位置的光源为焦点来消除阴影，垂直移动安装制动装置，不受手术房间的空间限制，使用安全可靠。

Adopts light sources from different positions for focusing to eliminate shadow from different parts of medical workers. Vertical mobile setup with braking device is not limited by the space of surgical room and is safe and reliable in use.

照度Lx	Illumination Lx	160000≥Ec≥50000
色温K	Color temperature	4000±500
显色指数	Color rendering index	100≥Ra≥85
光束深度	Depth of light beam(mm)	≥500

总辐照度值 (工作区域温升≤10℃) The value of the total irradiance (service range for temperature rise≤10 degree)	<6mW/m² lx
---	------------

总辐照度 (术者头部温升≤2℃) Total irradiance(Temperature rise on the head of patient under surgical operation≤2 degree):	≤600W/m²
---	----------

功率 Power	110V
灯泡额定电压 Rated voltage of bulb	AC 24V
灯泡额定功率 Rated power of bulb	125W
输入功率 Input power	150VA



AG-LT016

LED检查灯
LED Lamp



AG-LT016B

LED检查灯
LED Lamp



AG-LT019A-1

LED手术灯
LED Surgery Light



AG-AM001麻醉机
Anesthesia Machine



AG-AP001吸引器
Vacuum Suction Machine



AG-D0022高频X射线机(移动X光机)
High Frequency Mobile X-ray machine



AG-HXJ01呼吸机
Portable respirator



AG-DE001除颤仪
Defibrillation apparatus



AG-BZ009 十二道心电图机
12 Lead ECG Machine



AG-D0027洗片机
X-RAY Automatic Processor



AG-D0031吸引器
Vacuum Suction Machine



AG-KL001观片灯
X-ray film lamp



AG-D0033便携式呼吸机
Portable respirator



AG-D0037离心机
Lab centrifuge



AG-BZ013胎心仪
多普勒Fetal doppler



AG-900B注射泵
syringe pump



AG-DRS001高频数字X射线成像系统
High Frequency Digital Radiography System



AG-C101A04

LDR BED(Intelligent)

技术参数:

床体尺寸: 长度2350mm 宽度980mm

床面尺寸: 长度1990mm 宽度870mm

病床高度: 最低≤490mm 最高≥890mm (不含床垫)

背板上折角度: ≥65°

座板上折角度: ≥15°

床面后倾角度: ≥8°

脚板升降距离: 170mm

脚板外摆角度: ≥90°

脚板上折角度: ≥90°

护栏升降距离: 405mm

Bed dimension including bed end and bumper wheel: Length2350mm Width980mm

Bed dimension: Length1990mm Width870mm

Bed height: Lowest ≤490mm Highest≥890mm (exclude the cushion)

Back section turning: ≥65°

Seat section turning: ≥15°

Trendelenburg: ≥8°

Leg section lifting distance:170mm

Footrest swing out: ≥90°

Footrest folding: ≥90°

Guard rail up/down movement:405mm

Castor dia.:150mm

标准配件 Standard Accessories:

床垫 Bed cushion	1套(set)
污秽盆 Waste basin	1个(piece)
静脉输液架 IV Pole	1件(piece)
电源线 Power cord	1件(piece)

AG-C101A02B

手动产病床(普通) Manual LDR BED(Ordinary)

技术参数:

床尺寸包括床、保险杠轮:长2050毫米,宽980毫米

床尺寸:长1815毫米,宽880毫米

高度:580毫米~950毫米(排除垫)

背板角度调节: ≥70°

后倾: 5° 前倾: 12°

脚轮直径:125毫米

Technical Parameters:

Bed dimension including bed end and bumper wheel: Length 2050mm, width 980mm

Bed dimension: Length 1815mm, width 880mm

Min. & Max. table height: 580mm~950mm(exclude the cushion)

Back section turning: ≥70°

Trendelenburg and reverse trendelenburg: reverse trend≥5° trend≥12°

Leg section telescope: 550mm

Castor Diameter: 125mm



标准配件 Standard Accessories:

腿架 Leg holder	1对(pair)
把手 Grab handle	1对(pair)
脚踏板 Foot rest	1对(pair)
污秽盆 Waste Basin	1个(piece)
防水垫 Water proof cushion	1套(set)
输液架 I.V. Stand	1个(piece)
新生儿软垫 Cushion for newborn baby	1片(piece)
床垫 Cushion	1套 包括靠背垫、座垫、协助平台垫 (set Including back cushion, seat cushion and assist platform cushion)





AG-S108

妇科检查床
Gynecology exam table

床体尺寸(不包含拉手): 1400×720×880
背板上开角度: 0-70°
脚踏凳尺寸: 500mm长, 220mm高
脚踏凳行程: 450mm
辅助台行程: 420mm
脚蹬行程: 305mm
脚蹬外翻角度: 15°、30°、45°
最大病人承重: 240kg

Bed size(without handrails): 1400*720*880mm
Backrest upward angle: 0-70°
Foot-stool size: 500mm in length, 220mm in height
Foot-stool stroke: 450mm
Auxiliary table stroke: 420mm
Pedal stroke: 305mm
Foot pedal swing angle: 15°, 30°, 45°
Maximum loading capacity: 240kg



AG-C101A03

多功能电动产床 豪华妇科产床
超低位豪华产床
LDR BED(Low Starting Position)

技术参数:

床尺寸包括床,保险杠轮:长1930毫米,宽940毫米
床尺寸:长1815毫米,宽780毫米
高度:580毫米-860毫米(排除垫)

Bed dimension including bed end and bumper wheel: Length 1930mm, width 940mm

Bed dimension: Length 1815mm, width 780mm

Min. & Max. table height: 580mm-860mm(exclude the cushion)

标准配件 Standard Accessories:

搁腿架 Leg holder 1对 (pair)
污秽盆 Waste basin 1个 (piece)
输液架 I.V. Stand 1个 (piece)
手控 Hand control 1块 (piece)
电源线 Power cord 1根 (piece)

Back section turning: ≥60°
Seat section turning: ≥10°
Trendelenburg: ≥8°
Leg section up/down movement: 100mm
Footrest swing out: ≥45°
Footrest folding: ≥45°
Guard rail up/down movement: 360mm
Castor Diameter: 150mm



AG-C309

电动妇科椅
Electric Obstetric Bed

主要技术参数:

长: 2000mm
宽: 600mm
高: 550mm
升降行程 ≥400mm
背板上折 ≥70°
背板下折 ≥10°
臀板上折 ≥28°
臀板下折 ≥10°

Main specification:

Length:2000mm
Width:600mm
Height: 550mm
Adjustment: ≥400mm
Back section raised from the horizontal: ≥70°
Back section lowered from the horizontal: ≥10°
Seat section raised from the horizontal: ≥28°
Seat section lowered from the horizontal: ≥10°

配置清单:

主床1
手持操纵器1
麻醉屏架1
拉手2(带护套)
托腿架2(带垫)
搁臂板2(带垫)
污水盆2(不锈钢)

Standard Accessory:

Main bed 1
Hand control 1
Anesthesia frame 1
Handle 2 (with protecting bush)
Leg holder 2 (with cushion)
Arm holder 2 (with cushion)
Filtth basin 2 (stainless steel)





AG-S104A

妇科椅 (电动液压)
Gynecology Chair(Electric-Hydraulic)

技术参数:
台面尺寸:长1145毫米,宽500毫米
高度:550毫米~970毫米
背板角度调节:8°~50°
座位角度调节:≥25°
扶手外摆:0~90°
脚踏板外摆:≥45°
灯调整角度:60°
电压:50 hz交流220 v
标准配件
扶手1对
腿架1对
脚踏板1对
脚踏开关1件
无影灯一套

主部件为铝
附件和导轨为不锈钢

Technical Parameters
Table top dimension: length 1145mm, width 500mm
Min. & Max. table height: 550mm ~ 970mm
Back section turning: down ≥ 8° up ≥ 50°
Seat section turning: ≥ 25°
Arm rest swing out: 0 ~ 90°
Foot rest swing out: ≥ 45°
Lamp adjust angle: 60°
Power: AC 220V 50HZ
Standard accessories
Arm rest 1 pair
Leg holder 1 pair
Foot rest 1 pair
Foot switch 1 piece
Shadowless lamp 1 set



AG-S105C

多功能电动产床
Electric Gynecology Chair

床的尺寸/表面 Dimension of the bed /surface:
总长度 Overall length: 1250 ± 50mm
宽度 Width: 520 ± 20mm
高度 Height:
最高面高度 Max height of table top: 1200mm
最低面高度 Min height of table top: 480mm
座位板高度 Height of seating plate: 520mm

调整角度范围 Adjusting Scope of Tabletop:
前倾 Reversed trendelenburg: ≥ 20°
后倾 Trendelenburg: ≥ 5°



AG-S107

多功能电动产床
Electric Gynecology Chair

高: 480mm
升降行程 ≥ 590mm
背板上折 5° ~ 75°
臀板翻转角度 0~180°
背臀联动
腿架开关角度 0° ~ 30°

标配Standard accessories:
主床 Main Bed 1pc
脚踏开关 Foot switch 1pc

Height :480mm
Height adjustment ≥ 590mm
Backrest adjustment: 5° ~ 75°
Seat board turning: 0~180°
Seat Leg section on-off: 0° ~ 30°



AG-ECC01

手动检查床
Examination Couch (Back-rest Adjustable)
尺寸 Size: 1810 × 620 × 710mm (L × W × H)



AG-ECCZ1

电动检查床
Electric Examination Couch
尺寸 Size: 1900 × 700 × 480~700mm (L × W × H)



AG-ECC05

检查床
Examination Couch
尺寸 Size: 1800 × 600 × 700mm (L × W × H)



AG-ECC15

检查床
Examination Couch
尺寸 Size: 1800 × 600 × 700mm (L × W × H)
最小起订量 MOQ: 5 units



AG-ECC25

电动升降检查床
Electric Lift Massage Table
尺寸 Size: 1950 × 660 × 480~900mm (L × W × H)
最小起订量 MOQ: 5 units



AG-ECC24

液压检查床
Hydraulic Examination Couch
尺寸 Size: 1850 × 520 × 460~765mm (L × W × H)
最小起订量 MOQ: 5 units




AG-XD206

多功能电动医疗椅
Electric Medical Chair

技术参数 Technical data

配输液架, 餐桌 With IV stand, dinning table
 高度可调 Height can be adjustable
 尺寸 Size: L1880 x W920 x H1400 - 1640 mm
 背板角度调节 Back section: -10° - -80°
 腿板角度调节 Leg section: 10° - -80°
 座位角度调节 Seat segment angle: 0° - -30°
 前倾和后倾 Trendelenburg/Reverse trendelenburg: 12°
 最大承重量 Max loading weight :240kgs
 三个LINAK电机 Three linak motors
 选配: CPR,称重系统 Accessory: CPR, Weighting system



AG-XD107

Electric Blood Donation Chair

产品参数 Specification:

椅面尺寸 Chair dimension: 1700mm(曲线尺寸 Curve length) x 560mm
 椅面高度 Height: 570mm(min.)-830mm(max)
 头板折转角度 Head section turn: -15° ~ 25°
 背板折转角度 Back section lift: -5 ~ 70°
 腿板折转角度 Leg section turn: -65 ~ 5°
 踏板伸缩距离 Foot support telescopic: 100mm
 踏板折转角度 Foot support turn: 0-90°
 座板固定角度 Fixed angle of seat section: 5°

搁手板外摆角度 Arm board swing out: 0 ~ 45°
 搁手板上下调节距离 Arm board height adjustment: 100mm
 搁手板前后移动距离 Arm board slide: 100mm
 搁手板上下倾斜调节角度 Arm board declining angle adjustment: -15° ~ 10°
 搁手板面尺寸 Arm board dimension: 450 x 160mm
 前后倾斜角度 Trendelenburg and reverse trendelenburg: 0 ~ 12°
 扶手凹形弧度 Arm board concave radian: R108



AG-XD205

多功能手动医疗椅
Manual Blood Donation Chair

Seat width: 600mm
Seat width, including armrests: 880mm
Total length: 1880mm
Total weight: 70kg
Height: 530mm
Caster: dia.100mm
Back cushion: 190*490*120mm
Cushion: PU foam
Upholstery: PVC leather
Frame: Q235 steel
Safe maximum load: 240kg
Back segment angle: -12° -75°
Trendelenburg: -12°
Leg segment angle: -80° -12°
Color Optional: Green ,blue



AG-XD208B

多功能电动医疗椅
Electric Medical Chair

技术参数

尺寸:1880Lx600Wx530Hmm
背部调节: 0° ~75°
腿部调节: -80° ~0°
脚轮直径: 4x φ 100mm
面料: PU
重量: 85kg
最大承重: 200kg
2个LINK电机

Technical data

1. Size:1880Lx600Wx530Hmm
2. Back Section: 0° ~75°
3. Leg Section: -80° ~0°
4. Size of castors: 4x φ 100mm
5. Cushion: PU
6. Weight: 85kg
7. Max weight: 200kg
8. Two LINK Motors



AG-IIR001C UP-grade

高档婴儿培养箱
Infant Incubator

标准配置:

主体 (包括透明罩, 控制系统, 婴儿床, 水槽) 输液架, 储物柜, 静音脚轮

Standard configuration:

Main body(including the transparent hood, control system, infant bed, water trough), I.V.pole, mattress, cabinet, silence castors

可选功能

1. X射线盒式盘
2. 婴儿体重系统
3. 电动高度可调
4. LED光疗单元
5. 空气和氧气混合器
6. 双墙罩

Optional parts:

1.X ray cassette Tray
2.baby weight system
3.Electric height adjustable
4.LED phototherapy unit
5.Air and oxygen blender
6.Double wall hood

1.基于微处理器的伺服控温系统
2.LED温度显示屏
3.报警功能: 电源故障; 温度过高; 温度偏差; 温度传感器故障; 风扇故障
4. > 37℃温度设定
5.键盘锁、防止无意中更改设置
6.导热关闭功能
7.两档位湿度调节
8.双电源开关
9.婴儿床倾斜角度可调
10.推拉式水槽的清洗和加水
11.低噪音风扇, 使用寿命长
12.温度偏差和温度调节可在前面板上直接设置
13.静音报警按钮/重置按钮报警
14.六个操作窗口, 两侧的操作窗可旋转
15.孵化器内部的空气温度超过40℃时加热器自动断电
16.带四个抽屉柜

1. Microprocessor based servo controlled temperature system
2. Control modes: Air mode and baby skin mode
3. Alarms functions : Power failure; over temperature; temperature deviation; temperature sensor failure; fan failure
4. >37° C temperature setting
5. A keypad lock prevent inadvertent changes the setting
6. Humidity is adjustable continuously
7. Heating power displayed on panel
8. A second thermal cut-off function for more safety
9. The infant bed inclination angle is adjustable
10. Trough for pulling and pushing ,which is convenient to add the water and wash
11. Fan with low noise and long life-span
12. Temperature deviation and over temperature could be amended in front panel directly
13. Silence/reset button for alarms
14. six operating windows, two side operating windows can be swiveled
15. Heater is electrically cut-out if air temperature inside incubator exceed 40°C
16. Cabinet with four drawers

AG-IIR001A

婴儿暖箱
Infant Incubator

主要介绍:

1.微处理器控制的温度调节系统
2.LED温度显示屏
3.报警功能: 电源故障; 温度过高; 温度偏差; 温度传感器故障; 风扇故障
4. > 37℃温度设定
5.键盘锁、防止无意中更改设置
6.导热关闭功能
7.两个等级的湿度调节
8.板面显示加热功率
9.倾斜角度可调
10.推拉式水槽的清洗和加水

11.低噪音风扇, 使用寿命长
12.温度偏差和温度调节可在前面板上直接设置
13.静音报警按钮/重置按钮报警
14.六个操作窗口, 两侧的操作窗可旋转
15.带支架
16.当内部温度超过40℃时, 加热电源自动切断

Function features:

1.Temperature displayed by LED
2. Alarm functions: Power failure; Over temperature; Temperature deviation;
3.Temperature sensor failure; Fan failure;
4 >37 °C temperature setting
5 A keypad lock prevent inadvertent changes the setting
6 A second thermal cut-out function for more safety
7 Humidity by natural airway
8 Heating power displayed on panel
9 The infant bed inclination angle is adjustable





AG-IRW002B

高档婴儿保温台
Top Grade Infant Warmer

标准配置:

主体 (包括辐射源, 控制系统, 婴儿床, 支架) X线暗盒, 输液架
皮肤温度传感器, 托盘, 床垫, 透明保护板, 脚轮

Standard configuration

Main body (including the radiant source, control system, infant bed, bracket), X-ray cassette, I.V. Pole, skin temperature sensor, tray, mattress, transparent protector, castors.

- 1.基于微处理器的伺服控温系统
- 2.微处理器控制单元允许三个设置: 预热, 手动和自动
- 3.报警功能: 电源故障; 温度偏差; 温度传感器故障; 温度过高
- 4.APGAR计时器, 记录治疗时间
- 5.X射线暗盒
- 6.加热头可水平调节90°
- 7.婴儿床倾斜角度可调
- 8.可折叠透明护栏
- 9.采用防爆微晶石英红外辐射管为它的辐射源, 使加热迅速并且寿命长
- 10.可视听报警
- 11.第二热切断器功能更加安全
- 12.温度偏差可在面板上进行修正

Specification:

1. Three modes: pre-warm, manual and automatic
2. The microprocessor-controlled unit allows for three settings: Pre-warm, manual or automatic.
3. Alarm functions for safety includes: Power failure, temperature deviation, Temperature sensor failure, over temperature.
4. APGAR timer to record therapy time
5. Convenient X-ray cassette under infant mattress
6. Heater head can be adjusted ± 90 degree horizontally
7. Inclination of infant bed is adjustable
8. Transparent protector folded
9. Adopts anti-explosion micro crystal quartz infrared radiation tube for its radiation source, the heating is quick and life-span is long
10. Alarms are audible and visual
11. A second thermal cut-out function for more safety
12. Skin temperature sensor failure protect function avoid over temperature
13. Temperature deviation can be amended in front panel directly

AG-IRW003

新生儿保温台
Infant Radiant Warmer

- 1.基于微处理器的伺服控温系统
- 2.微处理器控制单元允许三个设置: 预热, 手动和自动
- 3.报警功能: 电源故障; 温度偏差; 温度传感器故障; 温度过高
- 4.APGAR计时器, 记录治疗时间
- 5.X射线暗盒
- 6.加热头可水平调节90°
- 7.婴儿床倾斜角度可调
- 8.观察灯角度可调
- 9.可折叠透明护栏
- 10.采用防爆微晶石英红外辐射管为它的辐射源, 使加热迅速并且寿命长
- 11.可视听报警
- 12.第二热切断器功能更加安全
- 13.温度偏差可在面板上进行修正
- 14.低压吸引器、氧气供应系统和LED光疗
- 15.电动调节功能

Brief introduction:

- 1.Microprocessor based servo controlled temperature system
- 2.Three modes: pre-warm,manual or automatic
- 3.Set temperature and skin temperature display separately
- 4.Special design inclined heater head is more convenience for X rays
- 5.Anti-explosion micro crystal quartz for radiation source
- 6.Self-testing function, various failure alarms by audible and visual
- 7.APGAR timer to record therapy time
- 8.Convenient X-ray cassette tray under infant mattress
- 9.Inclination of infant bed is adjustable
- 10.LED observe light angle is adjustable
- 11.Adopts anti-explosion micro crystal quartz infrared radiation tube for its radiation source, the heating is quick and life-span is long
- 12.Alarm functions for safety includes:Power failure, Temperature deviation, Temperature sensor failure,Over temperature.



AG-XB-Y1000
单管道输液泵
Single-channel Infusion Pump



AG-XB-Y1200
双管道输液泵
Double channel Infusion Pump



AG-BZ007
监护仪
Patient Monitor



AG-BZ014
多参数病人监护仪
Multi-parameter Patient Monitor



AG-BZ011

胎儿孕妇监护仪
Fetal monitor

Dimension: 360mm*310mm*120mm
12.1" color LED (0-90° rotation)
Resolution:1024*768 pixels
Type: Rechargeable lithium battery 14.8V/2200mAh
Charge Cycle: ≥500 times



AG-BZ012

胎儿孕妇监护仪
Fetal monitor

Dimension: 320mm*260mm*80mm
10.4" color LED (0-135° rotation)
Resolution:800*600 pixels
Type: Rechargeable lithium battery 7.4V/2200mAh
Charge Cycle: ≥500 times

AG-BU009

B超机
Laptop ultrasound scanner

Feature:

1. 15-inch LED screen
2. Windows embedded operating system
3. Built-in reporting workstation template (install the printer driver to print reports)
4. Built-in 8G SSD card hard disk
5. Built-in 8000 mA lithium battery (can work 3.5-4 hours)
6. Support each department inspection and measurement (routine examination: abdomen, obstetrics, gynecology, urology, small parts)
7. Reconstruction of three-dimensional features
8. Position marks: 140 kinds
9. With 44 kinds of pseudo-color display
10. Software upgrade: the software can be upgraded through U disk
- 11 Computer PC platform, built-in reporting workstation, micro-color, 3D reconstruction function



AG-BU005

多普勒扫描仪
Color Doppler Ultrasound System

Feature:

1. 15" LED Screen
2. High Definition image
3. 3D/4D imaging(option)
4. B,M,CF,PWD,PD,PDP,THI,TDI,real-time M,CWD
5. Abundant functions,Multilingual
6. Hardisk,USB storage
7. Color scanning rate:Screen degree: 85, full color vision
8. Color enhancements function:color Doppler energy,direction speed chart
9. Advanced and convenient direct network real-time image exchange,storage and display
10. Blood flow function





AG-WT006

ABS医疗工作车 ABS Medical Workstation Trolley

Height adjustment range for LCD and worksurface: 400 ± 20 mm, Lifting load ≤ 15 kg
 Worksurface height: 890mm-1290mm
 Large, smooth ABS worksurface dimensions:533mm(L) x455mm (W),Proposal Table load ≤ 5 kg
 Drawer inner dimensions: 379mm(L) x 298.5mm(W)x 84mm(H)
 Keyboard tray dimensions: 440mm(L) x 256mm(W)
 Mouse tray dimension: 208mm(L) x 200mm(W)
 Storage Bin with handle
 LCD column dimensions: 386mm(L)x80mm(W)x 64.5mm(H), Height adjustable with the help of internal lift construction.
 LCD hanger:height adjustment for LCD pivot is 13cm; 20° tilting for LCD; Applied for ≤ 22 inch LCD display.
 Proposal load ≤ 10 kg
 Base: PA6 + ABS cover
 Battery and UPS system inside the base; White LEDs provide light over keyboard tray
 Front-locking silent casters



AG-WT007

笔记本电脑车 ABS Medical Workstation Trolley

Notebook computer trolley
 Height adjustment range for worksurface: 400 ± 20 mm,
 worksurface height:750mm-1150mm,lifting load ≤ 15 Kg
 Large,smooth ABS worksurface dimension:533mm(L) x 455mm(W)
 Laptop cabinet dimension:430mm(L) x 346mm(W) x 56.5mm(H)
 Applied for ≤ 15 inch laptop
 Optional rotary keyboard tray, keyboard tray,dimension:440mm(L) x 256mm(W)
 Mouse tray dimension:208mm(L) x 200mm(W)
 Base:PA6 + ABS cover
 Silent medical casters



AG-MT014A

多功能电脑推车 Computer Trolley

技术参数:

1. 整车采用铝塑材质
2. 面板, 抽屉, 车身及托盘等由ABS工程塑料模具注塑成型
3. 左右两侧的立柱为铝合金材质
4. 整车台面配多功能护栏, 车的两侧配置物盒, 锐器盒, 金属挂篮, 开启式垃圾桶
5. 两大抽一小抽, 抽屉内配活动插片, 可自由组合
6. 台面以下配笔记本支架(电脑一体机及可充电电池)方便使用医院的管理系统
7. 及时查询病人信息
8. 中控锁, 静音脚轮(两个带刹车)

Technical Parameters

1. Made of high-strength ABS material
2. One-piece ABS plastic top board with arc handle and raised-edge design, covered transparent soft plastic glass
3. With four plastic-steel columns.
4. One small, two big drawers,inner with dividers can be organized easily and freely
5. Dust basket
6. Storage box
7. Disposable holder
8. Four luxurious noiseless casters, two with brakes

Package Size: 780 x 625 x 1170mm

Net Weight: 80kg
 Size: 660x520x1440mm



AG-WT002A

多功能电脑推车 Computer Trolley

技术参数:

1. 铝合金工作台面, 品牌电脑主机及升降系统
2. 显示屏, 品牌键盘鼠标, 高强度无线接收
3. ABS材质主机箱, 抽拉键盘托, ABS材质电池仓
4. 一体化铝合金底盘, 高强度, 耐碰撞
5. 配高强度静音脚轮, 带刹车

Technical Parameter

1. Tabletop is made of ABS material by use of one-time-molding technique
2. Fan-shaped design can add working area
3. Height of tabletop adjustable
4. Host Computer is built in, low power loss, D.C. stabilized source, include network
5. Battery is 12V-18.5V, high power and safe
6. With the space for mouse, scanner, keyboard
7. Fan-shaped chassis walk steadier
8. Four luxurious noiseless casters, two with brakes

Package Size: 560 x 600 x 1100mm

Net Weight:75kg

Size: 500 x 520 x 1330-1500mm

Function: Height adjust: 900-1070mm



DARK GREEN SERIES


AEGEAN产品设计师巧妙使用国际标准潘通色 (pantone) GG4004系列墨绿色编号: 7734UP
Designer of AEGEAN clever use Dark Green which is 7734UP in the international standard color of PANTONE GG4004 series.

设计主题—给医院色彩
AEGEAN产品设计理念 “简约而不简单”
Design Theme—Give life and vigor to hospital
AEGEAN design philosophy--- ‘simple but not easy’



AG-ET001A1

医用急救推车
Emergency Trolley
AG-ET001A1(850mm) 尺寸 Size: 850 × 520 × 950mm
AG-ET001A1(750mm) 尺寸 Size: 750 × 480 × 920mm



AG-GS001
尺寸 Size:940 × 620 × 450mm



AG-GS002
尺寸 Size:940 × 620 × 460mm



AG-GS003
尺寸 Size:940 × 620 × 460mm



AG-GS004
尺寸 Size:950 × 630 × 435mm



AG-GS005
尺寸 Size:970 × 640 × 460mm



AG-GS006
尺寸 Size:970 × 660 × 450mm



AG-GS007
尺寸 Size:910 × 510 × 435mm



AG-GS008
尺寸 Size:950 × 645 × 470mm



AG-GS009
尺寸 Size:950 × 645 × 470mm



AG-GS010
尺寸 Size:910 × 740 × 410mm

AG-MT011A1

双面药车
Medicine Trolley
尺寸 Size: 850 × 520 × 950mm





AG-MT006 医用药车
Medicine Trolley
尺寸 Size: 590 x 470 x 920mm



AG-MT001A1 医用药车
Medicine Trolley
尺寸 Size: 850 x 520 x 950mm



AG-ET007B3 急救推车
Emergency Trolley
尺寸 Size: 720 x 470 x 920mm



AG-AC001
陪护椅
Accompany chair
折叠尺寸 Folded size:
620 x 720 x 940mm
展开尺寸 Unfolded size:
620 x 1920mm



AG-AC012
豪华陪护椅
Accompany chair
折叠尺寸 Folded size:
830 x 850 x 420mm
展开尺寸 Unfolded size:
2100 x 850 x 610mm



AG-AC005
豪华陪护椅
Accompany chair
折叠尺寸 Folded size:
650 x 780 x 1020mm
展开尺寸 Unfolded size:
2060 x 650 x 450 mm



AG-AT001A3 麻醉车
Anesthesia Trolley
尺寸 Size: 750 x 480 x 1340mm



AG-NT003A1 治疗车
Nursing Trolley
尺寸 Size: 850 x 520 x 950mm



AG-LPT006B 治疗车
Clinical Trolley
尺寸 Size: 750 x 480 x 1000mm



AG-AC007
豪华陪护椅
Luxury Accompany Chair
Folded: 800 x 710 x 430 mm
(back section height 980mm)
Unfolded size: 2040 x 710 x 430 mm
Specification: PVC, Powder coating steel ,
solid wooden handrail



AG-TC001
豪华输液椅
Luxury Transfusion Chair
尺寸 size:
770 x 680 x 1080mm



AG-AC006
豪华输液椅
Luxury Transfusion Chair
尺寸 size:
730 x 840 x 1050mm



AG-CHT003 病历夹推车
Patient Record Trolley
尺寸 Size: 750 x 410 x 950mm



AG-MT032 治疗车
Treatment Trolley
尺寸 Size: 620 x 420 x 930mm



AG-UTB09 治疗车
Multi-function Cart
尺寸 Size: 640 x 445 x 1010mm



AG-AC009
陪护椅
Accompany chair
折叠尺寸 Folded size:
660 x 680 x 410mm
展开尺寸 Unfolded size:
1860 x 680 x 400mm



AG-AC011
沙发床
Sofa bed
Unfolded Size:
Length:1920mm,Height:260mm.
Folded Size:
Length:800mm x 1050mm x 790mm,seat height:380mm
Color: Lake blue
Material: Fine linen fabrics, sponges, steel



AG-TC003
输液椅
Infusion Chair
Size:
890 x 830 x 1050mm
Angle adjustable:
110 - 160°



AG-BC001 Bedside Locker
尺寸 Size: 480 × 480 × 750mm



AG-BC005 Bedside Locker
尺寸 Size: 475 × 475 × 750mm



AG-BC006C Bedside Locker
尺寸 Size: 480 × 480 × 760mm



AG-BC008 Bedside Locker
尺寸 Size: 510 × 500 × 700mm



AG-BC010 Bedside Locker
尺寸 Size: 460 × 440 × 780mm



AG-BC022 Bedside Locker
尺寸 Size: 500 × 500 × 780mm



AG-BC023 Bedside Locker
尺寸 Size: 500 × 500 × 780mm



AG-BC025 Bedside Locker
尺寸 Size: 600 × 520 × 700mm



AG-NS001 Doctor Stool
尺寸 Size: 400 × 380-510mm



AG-NS005B Doctor Stool
Seat Size: 430 × 400mm
Seatback Size: 400 × 300mm



AG-OBT002 Over Bed Table
尺寸 Size: 900 × 450 × 720-1020mm



AG-OBT009 Over Bed Table
尺寸 Size: 900 × 450 × 720-1020mm



AG-OBT012 Over Bed Table
尺寸 Size: 800 × 400 × 800-1110mm



AG-NS004 Patient Chair
尺寸 Size: 268 × 338 × 455mm



AG-SS003 S.S Instrument Cupboard
尺寸 Size: 900 × 400 × 1800mm



AG-SS004 Instrument Cupboard
尺寸 Size: 900 × 400 × 1800mm



AG-SS007 Fan-shape Instrument Trolley
尺寸 Size: 1370 × 400 × 1000mm



AG-SS008 S.S Mayo Cart
尺寸 Size: 660 × 400 × 940-1400mm



AG-KB001 Kick Bucket
尺寸 Size: ϕ 23/ ϕ 29 × 24CM



AG-SS013 S.S Dirty Linen Trolley
尺寸 Size: 700 × 500 × 900mm



AG-SS020 Instrument Trolley
尺寸 Size: 850 × 500 × 950mm



AG-SC001 Bed screen
尺寸 Size: L2000 × H1800mm × W500



AG-SS019B SS Linen Trolley
尺寸 Size: 780 × 570 × 880mm



AG-WAS009 Water Sinks
尺寸 Size: 1300*630 / 450*1800mm



AG-TWC002 Waiting chair
尺寸 Size: 1750 × 680 × 770mm



AG-TWC004 Waiting chair
尺寸 Size: 1910 × 750 × 800mm



AG-SS010A Linen Trolley
尺寸 Size: 900 × 500 × 1050mm



AG-SS032 S.S Flat Cart
尺寸 Size: 850 × 500 × 900mm



AG-SS034 S.S Dressing Cart
尺寸 Size: 380 × 450 × 920mm



AG-SS042 S.S Instrument Trolley
尺寸 Size: 850 × 500 × 900mm



AG-FS008 Double foot steps
尺寸 Size: 410 × 230 × 280mm



AG-SS080 S.S Stack Instrument Trolley
尺寸 Size: 1200 × 600 × 1050mm



AG-SS042D Treatment Trolley
尺寸 Size: 690 × 510 × 880mm



AG-SS066 S.S Oxygen Cylinder Cart
尺寸 Size: 330 × 330 × 1200mm