



We are near grasslands, forest, ocean...

TÜVRheinland® CE 0197
Precisely Right.



Spectrum For Medical And Laboratory Co.,Ltd

- Jordan - Amman - Khalda

www.spectrum-medlab.com
info@spectrum-medlab.com



Spectrum[®]

MEDICAL AND LAB SUPPLIES

OXYGEN THERAPY

Oxygen therapy is a treatment that provides you with extra oxygen, a gas that your body needs to work well. Normally, your lungs absorb oxygen from the air. However, some diseases and conditions can prevent you from getting enough oxygen. Oxygen therapy may help you function better and be more active. Oxygen is supplied in a metal cylinder or other container. It flows through a tube and is delivered to your lungs in one of

the following ways:

1. Through a nasal cannula, which consists of two small plastic tubes, or prongs, that are placed in both nostrils.
2. Through a face mask, which fits over your nose and mouth.
3. Through a small tube inserted into your windpipe through the front of your neck.

All materials we used in construction of the oxygen mask, and the oxygen tubing are soft and smooth surface without sharp edge and object. They have no undesirable effects on the oxygen/medication passing through under ordinary conditions of use. Mask material is hypoallergenic and shall resist ignition and rapid burn.



Spectrum for medical and lab supplies is an innovative sourcing & distribution company, involved in exclusive representation and global supply of a wide range of medical devices, commodities, and high-quality pharmaceuticals

Our leading products are Respiratory, Anesthetic Medical Products, Laboratory Products, etc. The company has passed ISO13485, TUV CE certification of more than 40 products, and a number of patents.

The founders of this company and their staff have worked for the Iraqi ministry of health and occupied high rank positions in the government whether that was for a short period or long it results in know how when it comes to successfully winning and fulfilling tenders and contracts

We will continue regarding laws and regulations as cornerstone, implementing full staff management, pursuing technological innovation and excellent quality.

Under the concept of "people first", our company has formed an excellent team with great team spirit, full of vigor and great execution ability. We warmly welcome customers both at home and abroad to consult and cooperate. We are sure to provide best service for every customer sincerely!



LB1110 SIMPLE OXYGEN MASK



Made of clear, medical grade PVC
with adjustable nose clip, elastic strap with/without latex
Medium concentration, approx 50% .
Color: green transparent, white transparent and light blue transparent is available.
Length: with 2m tubing, different length can be customized,
Anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Cat.No.	Size
LB111000	XS(Neonate)
LB111001	S(Pediatric Standard)
LB111002	M(Pediatric Elongated)
LB111003	L(Adult Asian)
LB111004	L+(Adult European)
LB111005	XL(Adult Elongated)

LB1210 VENTURI MASK



Made of clear, medical grade PVC
With adjustable nose clip, elastic strap with/without latex
Two color coded diluters for seven concentrations
Green diluter for 24%, 26%, 28%, 30%, White diluter for 35%, 40%, 50%
Color: green transparent, white transparent and light blue transparent is available.
Length: with 2m tubing, different length can be customized,
Anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Cat.No.	Size
LB121000	XS(Neonate)
LB121001	S(Pediatric Standard)
LB121002	M(Pediatric Elongated)
LB121003	L(Adult Asian)
LB121004	L+(Adult European)
LB121005	XL(Adult Elongated)

Table of correspondent oxygen concentration for different oxygen flow

Adjusted Oxygen Concentration	Oxygen Flow
24%	3L/min
26%	3L/min
28%	6L/min
30%	6L/min
35%	9L/min
40%	12L/min
50%	15L/min

LB1220 ADJUSTABLE VENTURI MASK

Made of clear, medical grade PVC
 With adjustable nose clip, elastic strap with/without latex
 Single diluter for one concentration, six diluters can be selected
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized,
 Anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Cat.No. Type A	Cat.No. Type B	Size
LB122000A	LB122000B	XS(Neonate)
LB122001A	LB122001B	S(Pediatric Standard)
LB122002A	LB122002B	M(Pediatric Elongated)
LB122003A	LB122003B	L(Adult Asian)
LB122004A	LB122004B	L+(Adult European)
LB122005A	LB122005B	XL(Adult Elongated)




Table of correspondent oxygen concentration for different oxygen flow (two types)

LB1220A: A type with below 6pcs oxygen concentrations

Adjusted Oxygen Concentration	Oxygen flow
24%(blue)	4L/min
28%(yellow)	4L/min
31%(white)	6L/min
35%(dark green)	8L/min
40%(pink)	8L/min
50%(orange)	10L/min

LB1220B: B type with below 7pcs oxygen concentrations

Adjusted Oxygen Concentration	Oxygen flow
24%(blue)	4L/min
28%(yellow)	4L/min
31%(white)	6L/min
35%(dark green)	8L/min
40%(pink)	8L/min
50%(orange)	10L/min
60%(red)	15L/min

LB1310 OXYGEN MASK WITH RESERVOIR BAG

Made of clear, medical grade PVC
 With reservoir bag
 With adjustable nose clip, elastic strap with/without latex
 Oxygen mask with or without check valve
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized,
 Anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.



Cat.No.	Size
LB131000	XS(Neonate)
LB131001	S(Pediatric Standard)
LB131002	M(Pediatric Elongated)
LB131003	L(Adult Asian)
LB131004	L+(Adult European)
LB131005	XL(Adult Elongated)

LB1510 TRACHEOSTOMY MASK

Made of clear, medical grade PVC
 With elastic strap with/without latex
 Swivel 360° adaptor with 18f/22m port,
 Excellent compatible with corrugated tubing & other connection device.
 Packing: packed in individual PE pack. 50pcs/ctn or 100pcs/ctn

Cat.No.	Size
LB151001	M(Child)
LB151002	L(Adult)



LB1410 NASAL OXYGEN CANNULA



NASAL OXYGEN PRONG



LB1401
Curved Prong

LB1402
Straight Prong

LB1403
Soft Prong

Made of clear, medical grade PVC
 Soft prong, standard prong, curved prong is available
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized, anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Size	Cat.No. : LB1411 (With Curved Prong)	Cat.No.: LB1412 (With Straight Prong)	Cat.No. : LB1413 (With Soft Prong)
XXS(Neonate)	/	LB141201	LB141301
XS(Infant)	/	LB141202	LB141302
S(Child)	LB141103	LB141203	LB141303
L(Adult)	LB141104	LB141204	LB141304

LB1420

NASAL OXYGEN CANNULA WITH EAR GUARDS



Made of clear, medical grade PVC
 With medical grade silicone rubber earloop
 Soft prong, standard prong, curved prong is available
 Color: green transparent, white transparent
 and light blue transparent is available.
 Length: with 2m tubing, different length can be customized,
 Anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack.
 50 pcs/carton or 100 pcs/carton.

Size	Cat.No. : LB1421 (With Curved Prong)	Cat.No.: LB1422 (With Straight Prong)	Cat.No. : LB1423 (With Soft Prong)
XXS(Neonate)	/	LB142201	LB142301
XS(Infant)	/	LB142202	LB142302
S(Child)	LB142103	LB142203	LB142303
L(Adult)	LB142104	LB142204	LB142304

LB1120


CO₂ SAMPLING MASK



Made of clear, medical grade PVC
 With adjustable nose clip, elastic strap with/without latex
 Medium concentration
 The luer lock of CO₂ sampling line can be male or female.
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized,
 Anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.


Cat.No.	Size
LB112000	XS(Neonate)
LB112001	S(Pediatric Standard)
LB112002	M(Pediatric Elongated)
LB112003	L(Adult Asian)
LB112004	L+(Adult European)
LB112005	XL(Adult Elongated)

LB1430 ETCO₂ SAMPLING CANNULA



Just Adult Size is available with curved prong.
Designed for oxygen supply and sampling of exhaled CO₂ gas simultaneously.
By delivering oxygen through one prong and sampling exhaled gas from the other,
the nasal cannula can provide end tidal values comparable of those achieved with intubated patients.
The Connector for CO₂ sampling line can be male or female connector.
And the connector for O₂ supply can be funnel connector or luer lock universal connector

LB1610 OXYGEN CONNECTING TUBE



Made of clear, medical grade PVC
With connector
Color: green transparent, white transparent
And light blue transparent is available.
Length: with 2m tubing, different length can be customized,
Anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack.
50 pcs/carton or 100 pcs/carton.

LB1620 RESPIRATORY GAS MONITORING LINE



Respiratory gas monitoring lines in various lengths designed to fit between the patient end of the breathing system, via a suitable luer lock connection, and the anesthesia machine respiratory gas monitor connection.
Made of clear, medical grade PVC
Color: green transparent, blue transparent, white transparent are available.
The connector can be male or female.
Packing : packed in individual pack,
50pcs or 100pcs/carton.



AEROSOL THERAPY

Aerosol therapy is the delivery of medication during the inhalation phase of respiration.

Drugs are administered using aerosol therapy for the following reasons:

1. The ability to use a smaller amount of drug in a greater concentration.
2. The ability to deliver a drug to a specific site.
3. The ability to reduce the potential for systemic side effects.

Patients may use aerosol therapy for a few days or a lifetime, depending on their illness or disease state. A patient with an acute respiratory illness may use a single drug with a single device for a short period of time, whereas a patient with chronic respiratory illness such as asthma may use multiple drugs and multiple devices over an extended period of time.

Our Oxygen Tubing are latex free, soft and smooth surface without sharp edge and object, They are no undesirable effects on the Oxygen/Medication passing through under ordinary conditions of use. Mask Material are hypoallergenic and shall resist ignition and rapid burn.



LB20A NEBULIZER JAR

Volume: 6cc
Under the condition of gas flow: 7L/Min, atomizing volume is 200mg/min. In every 200mg, the percentage of particulate with the diameter less than 5um is more than 70%.

CODE	VOLUME	GAS FLOW (LPM)	TOR (mg/min)	Φ<5um (%)
LB20A	6CC	7	200	70%




LB20B NEBULIZER JAR

Volume: 6cc
Under the condition of gas flow: 7L/Min, atomizing volume is 280mg/min. In every 280mg, the percentage of particulate with the diameter less than 5um is more than 70%.

LB20B Patent Number : ZL 2014 2 0766313.3

CODE	VOLUME	GAS FLOW (LPM)	TOR (mg/min)	Φ<5um (%)
LB20B	6CC	7	280	70%



LB20C NEBULIZER JAR

Volume: 10cc
Under the condition of gas flow: 7L/Min, atomizing volume is 250mg/min. In every 250mg, the percentage of particulate with the diameter less than 5um is more than 70%.

LB20C Patent Number : ZL 201520637654.5

CODE	VOLUME	GAS FLOW (LPM)	TOR (mg/min)	Φ<5um (%)
LB20C	10CC	7	250	70%



LB20D NEBULIZER JAR (Adjustable)

Volume: 6cc
Under the condition of gas flow: 7L/Min, atomizing volume is 150-450mg/min. The percentage of particulate with the diameter less than 5um is more than 70%.

LB20D Patent Number : ZL 201520289256.9

CODE	VOLUME	GAS FLOW (LPM)	TOR (mg/min)	Φ<5um (%)
LB20D	6CC	7	150-450	70%

LB2011 NEBULIZER MASK



Made of clear, medical grade PVC
 Operates even at horizontal position capacity: 6cc
 With adjustable nose clip, elastic strap with/without latex
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized,
 Anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Size	Cat.No
XS(Neonate)	LB201100
S(Pediatric Standard)	LB201101
M(Pediatric Elongated)	LB201102
L(Adult Asian)	LB201103
L+ (Adult European)	LB201104
XL(Adult Elongated)	LB201105



LB20A



LB20B



LB20C

LB20D

* Different type of nebulizer jars are available.

LB2021 NEBULIZER MASK WITH 360° SWIVEL CONNECTOR



Made of clear, medical grade PVC
 Swivel 360° connector
 Operates even at horizontal position capacity: 6cc
 Without adjustable nose clip, elastic strap with/without latex
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized,
 Anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Size	Cat.No.
M(Child)	LB202102
L(Adult)	LB202103



LB20A



LB20B



LB20C



LB20D

* Different type of nebulizer jars are available.

LB2012 NEBULIZER KIT



LB20A

LB20B

LB20C

LB20D

* Different type of nebulizer jars are available.

Made of clear, medical grade PVC
 Operates even at horizontal position capacity: 6cc
 Contain: 1pc adult mask and 1pc child mask,
 T connector, mouthpiece, Nebulizer Jar and 2m tubing
 With adjustable nose clip, elastic strap with/without latex
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized,
 Anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

LB2013 NEBULIZER KIT WITH MOUTH PIECE AND CORRUGATED TUBE



LB20A

LB20B

LB20C

LB20D

* Different type of nebulizer jars are available.

Made of clear, medical grade PVC
 Operates even at horizontal position capacity: 6cc
 With T-piece and mouthpiece
 Color: green transparent, white transparent and light blue transparent is available.
 Length: with 2m tubing, different length can be customized, anti-crush tube can ensure oxygen follow even if the tube is kinked.
 Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

**LB2014
NEBULIZER KIT
WITH MOUTH PIECE**



LB20A

LB20B

LB20C

LB20D

* Different type of nebulizer jars are available.

Made of clear, medical grade PVC
Operates even at horizontal position capacity: 6cc
With T-piece and mouthpiece
Color: green transparent, white transparent and light blue transparent is available.
Length: with 2m tubing, different length can be customized, anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

**LB2015
NEBULIZER KIT
WITH MOUTH PIECE**



LB20A

LB20B

LB20C

LB20D

* Different type of nebulizer jars are available.

Made of clear, medical grade PVC
Operates even at horizontal position capacity: 6cc
Color: green transparent, white transparent and light blue transparent is available.
Length: with 2m tubing, different length can be customized,
Anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.



LB2511 NEBULIZER MASK

Made of clear, medical PVC
Size: M(Child), L (Adult)
Volume: 8cc
Packing: packed in individual pe pack



LB3810 Bite Block with Strap

Specifications

- To prevent biting gastroscopy channel
- For gastroscopy test and endoscopic operation



LB382002

LB382001

LB3820 Nose Clip

Specifications:

- Elasticized frame to press the soft tissue of nose wing
- 100% latex free
- Two material (TPE and POM) option
- Disposable, single use, non-sterile
- Color: White, Blue

Cat.No.	Specification	Color
LB382001	Material:TPE	White
LB382002	Material:Foam	Blue



LB 3830 Mouth Opener

Features:

The mouth opener is intended for opening mouth in the first-aid.
Material :LDPE

LB2120 SPACER FOR AEROSOL




Spacer Body Size : 50x153 (bottom diameter x height)
Volume : 175ml
Mouthpiece : 22.5x15.5mm(top of mouthpiece dimension)
Spacer Body Weight: 50g
Mask Weight : Small: 19.7G, Medium:23.7G, Large:43.8G
Material : 100% Latex free, Silicone Rubber, Anti-static PETG and TPR.

Cat.No.	Size	Applicable People	Color
LB212001	S	0-18 Months	Red
LB212002	M	1-5 Years	Yellow
LB212003	L	5+ Years	Purple/blue

AIRWAY MANAGEMENT

Airway management is a set of medical procedures performed in order to prevent airway obstruction and thus ensure an open pathway between a patient's lungs and the outside world. This is accomplished by clearing or preventing obstructions of airways, often referred to as choking, caused by the tongue, the airways themselves, foreign bodies or materials from the body itself, such as blood or stomach content, the latter resulting in aspiration. Airway management can be divided into two categories: Basic and advanced airway management. Basic techniques are simple to perform even by non-health care professionals and do not require use of medical equipment, whereas advanced techniques require special training and medical equipment. Advanced airway management is further categorized in increasing order of invasiveness into supraglottic devices, such as oropharyngeal or nasopharyngeal airways, followed by infraglottic techniques, such as tracheal intubation, and finally surgical methods. Airway management is a primary consideration in cardiopulmonary resuscitation, anaesthesia, emergency medicine, intensive care medicine and first aid.



LB3011 ANESTHESIA MASK

Made from Latex/DEHP free PVC material.
Air cushion assures comfortable face fitting.
Tail Check valve/Upright Check valve/Without check valve type are available.
Single Patient use.
Color-Coded hook options.

Cat.No.	Size	Color	Connector
LB301100	0# - Neonate	N/A	15mm OD
LB301101	1# - Infant	N/A	15mm OD
LB301102	2# - Pediatric	Yellow	22mm ID
LB301103	3# - Adult-Small	Green	22mm ID
LB301104	4# - Adult-Medium	Red	22mm ID
LB301105	5# - Adult-Large	Blue	22mm ID
LB301106	6# - Adult-Extra Large	Orange	22mm ID



LB3021 SOFT ANESTHESIA MASK

Made from Latex/DEHP free PVC material.
Soft adjustable face cushion assures comfortable face fitting.
With Upright Check valve.
Single Patient use.
No scents or with Vanilla scents are available.

Cat.No.	Size
LB3021000	Size00# - Premature
LB302100	size0# - Neonate
LB302101	size1# - Infant
LB302102	size2# - Pediatric
LB302103	size 3# - Adult-small
LB302104	size 4# - Adult-medium
LB302105	size 5# - Adult-large



LB3012 PVC FREE ANESTHESIA MASK

Made from TPE and PP Material, contain no PVC and are DEHP free
The product provides a comfortable superior airtight fit with a clear shell to enhance patient face visibility.
This product ancillary use with anesthesia breathing circuit.
Non-reusable, single use only

Cat.No.	Size	Color	Connector
LB301201	size1# - Neonate	Blue	15mm OD
LB301202	size2# - Infant	Orange	15mm OD
LB301203	size3# - Pediatric	Pink	22mm ID
LB301204	size4# - Adult-S	White	22mm ID
LB301205	size5# - Adult-M	Yellow	22mm ID
LB301206	size6# - Adult-L	Green	22mm ID



**LB3030
SILICONE ANESTHESIA MASK (ONE-PIECE)**

All-in-one piece soft silicone cushion for better seal and comfort.
Ergonomic shape for lower dead space and better air circulation.
Made of 100% clear, medical grade silicon, latex free
Reusable, all parts are fully autoclavable up to 134°C
Packing: Bulk or PE pack available

Cat.No.	Size	Connector
LB303000	0# - Neonate	15mm OD
LB303001	1# - Infant	15mm OD
LB303002	2# - Child	22mm ID
LB303003	3# - Adult-small	22mm ID
LB303004	4# - Adult-medium	22mm ID
LB303005	5# - Adult-large	22mm ID
LB303006	6# - Adult-extra large	22mm ID

**LB3031
SILICONE ANESTHESIA MASK (ROUNDED)**

All-in-one piece soft silicone cushion for better seal and comfort.
Ergonomic shape for lower dead space and better air circulation.
Made of 100% clear, medical grade silicon, latex free
Reusable, all parts are fully autoclavable up to 134°C
Packing: Bulk or PE pack available

Cat.No.	Size	Connector
LB303100	0# - Neonate	15mm OD
LB303101	1# - Infant	15mm OD
LB303102	2# - Child	22mm ID

**LB3032
SILICONE ANESTHESIA MASK (TWO-PIECE)**


Soft silicone cushion encased in hard shell for ease to handle and better seal.
Made of 100% clear, medical grade silicon and pc material, latex free.
This products ancillary use with disposable anesthesia breathing circuit.
Reusable, all parts are fully autoclavable up to 134°C
Packing: Bulk or PE pack available

Cat.No.	Size	Connector
LB303203	3# - Adult-small	22mm ID
LB303204	4# - Adult-medium	22mm ID
LB303205	5# - Adult-large	22mm ID



**LB3041
Silicone CPAP Mask and Headgear**

Material: medical silicone
Size: M, L
Sealing edge of contour line and head belt eliminate discomfort caused by the oppression to facial nerve
The design of nasal part fits well with the nose and leaves little air blind angle, which limits breath of exhaust
And increases comfort
durable silicone skirt, five points head belt, anti-fog design of standard nose and mouth cover



**LB 3042
CPAP Mask**

CPAP Mask is an interface for the treatment of obstructive sleep apnea using continuous positive airway pressure.
It Provide a good seal when fitted to the patient whilst providing maximum comfort.
The Accessories (Optional): PEEP Valve and Headgear

Cat.No.	Size	Connector
LB304101	Adult-S	22mm ID
LB304102	Adult-M	22mm ID
LB304103	Adult-L	22mm ID

LB3052 Accessories Headgear
Size: Pediatric /Adult

LB3051
Accessories
PEEP Valve



LB3110 SILICONE RESUSCITATOR

LB3110L LIQUID SILICONE RESUSCITATOR

LB3120 SEBS RESUSCITATOR

LB3130 PVC RESUSCITATOR



Silicone resuscitator is intended for pulmonary resuscitation
 The set includes oxygen reservoir bag, silicone mask and oxygen tubing
 Reusable, autoclave 134°C helping to prevent cross infection and contamination
 Guaranteed for 5 years of alternatively 20 times autoclave cycles
 Patient connector is 22/15mm.

SEBS/PVC resuscitator is disposable and intended for pulmonary resuscitator
 The set includes oxygen reservoir bag, PVC mask, and oxygen tubing
 100% medical level SEBS/PVC material, latex-free components
 All-in-one intake valve connects to oxygen reservoir bag directly.
 Patient connector is 22/15mm

Cat.No. Silicone	Cat.No. Liquid Silicone	Size	Specification	Package
LB311001	LB311001L	Adult: Adults and Children	Reusable adult 1650ml resuscitator(W/60cm H ₂ O POP-OFF) Silicon mask 5# 2000ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:60cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available
LB311002	LB311002L	Pediatric: Children and Infants	Reusable pediatric 600ml resuscitator(W/40cm H ₂ O POP-OFF) Silicone mask 2# 1600ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:40cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available
LB311003	LB311003L	Infant: Infants and Neonates	Reusable infant 280ml resuscitator(W/40cm H ₂ O POP-OFF) Silicone mask 1# for infant Or silicone mask 0# for neonate. 1600ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:40cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available

Cat.No. SEBS	Cat.No. PVC	Size	Specification	Package
LB312001	LB313001	Adult: Adults and Children	Single use adult 1650ml resuscitator(W/60cm H ₂ O POP-OFF) PVC anesthesia mask 4# 2000ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve: 60cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available
LB312002	LB313002	Pediatric: Children and Infants	Single use Pediatric 600ml resuscitator(W/40cm H ₂ O POP-OFF) PVC anesthesia mask 2# 1600ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:40cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available
LB312003	LB313003	Infant: Infants and Neonates	Single use infant 280ml resuscitator(W/40cm H ₂ O POP-OFF) PVC anesthesia mask 1# for infant Or PVC anesthesia Mask 0# for neonate 1600ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:40cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available



**LB3210
DISPOSABLE SILICONE LARYNGEAL MASK**

The tube is made from 100% medical silicone, No poison and harmless
Sterilization: EO
Packing: 1pc/blister packing

Cat.No.	Size	Applicable People
LB321010	1#	<5kg
LB321015	1.5#	5-10kg
LB321020	2#	10-20kg
LB321025	2.5#	20-30kg
LB321030	3#	30-50kg
LB321040	4#	50-70kg
LB321050	5#	>70kg



**LB3220
REUSABLE SILICONE LARYNGEAL MASK**

The tube is made from 100% medical silicone, No poison and harmless
Sterilization: EO
Packing: 1pc/blister packing

Cat.No.	Size	Applicable People
LB322010	1#	<5kg
LB322015	1.5#	5-10kg
LB322020	2#	10-20kg
LB322025	2.5#	20-30kg
LB322030	3#	30-50kg
LB322040	4#	50-70kg
LB322050	5#	>70kg



**LB3230
DISPOSABLE PVC LARYNGEAL MASK**

The tube is made from 100% medical PVC, Non-toxic PVC
No poison and harmless
Sterilization: EO
Packing: 1pc/blister Packing

Cat.No.	Size	Applicable People
LB323010	1#	<5kg
LB323015	1.5#	5-10kg
LB323020	2#	10-20kg
LB323025	2.5#	20-30kg
LB323030	3#	30-50kg
LB323040	4#	50-70kg
LB323050	5#	>70kg

**LB3310
GUEDEL AIRWAYS**

Made from Medical semirigid Non-toxic PE material.
Color-Coded bite block is designed for easy identification and prevent biting down so that can avoid blocking airway (PP material).
Guedel style enclosed center channel
Packing: Packed in individual PE bag or blister packing

Cat.No.	Code	Size(mm)	Color
LB331040	000#	40	Pink
LB331050	00#	50	Blue
LB331060	0#	60	Black
LB331070	1#	70	White
LB331080	2#	80	Green
LB331090	3#	90	Yellow
LB331100	4#	100	Red
LB331110	5#	110	Light Blue
LB331120	6#	120	Orange



**LB3320
BERMAN AIRWAYS**

Atraumatic one piece design for optimum patient comfort and safety.
Rigid bite block to ensure airway potency.
Two side channel design provides passage of suction catheter to clear secretions and also provides space for oral care.
Packing: Packed in individual PE bag or blister packing.

Cat.No.	Code	Size(mm)
LB332040	LB333040	000# 40
LB332050	LB333050	00# 50
LB332060	LB333060	0# 60
LB332070	LB333070	1# 70
LB332080	LB333080	2# 80
LB332090	LB333090	3# 90
LB332100	LB333100	4# 100
LB332110	LB333110	5# 110



LB3320



LB3330

**LB3410
NASOPHARYNGEAL AIRWAY**

Made from soft PVC.
Three types of connectors are available.
For nasal airway management only.

Cat.No.	Size ID(mm)
LB341025	2.5
LB341030	3.0
LB341035	3.5
LB341040	4.0
LB341045	4.5
LB341050	5.0
LB341055	5.5
LB341060	6.0
LB341065	6.5
LB341070	7.0
LB341075	7.5
LB341080	8.0
LB341085	8.5
LB341090	9.0
LB341095	9.5



CPR PRODUCTS

Combination of rescue breathing and chest compressions delivered to victims thought to be in cardiac arrest when cardiac arrest occurs, the heart stops pumping blood. CPR can support a small amount of blood flow to the heart and brain to "buy time" until normal heart function is restored.

LB3510 CPR MASK



Contents: foldable cushion mask with oxygen port, one-way valve
With pp box, box color has : orange, blue, green, yellow, white


LB3610 CPR BREATHING FILM



LB361001
With Blue Pig Valve

LB361002
With White Valve

LB362001
With Filter



Straight

Curved

Handle

LB3710 ANESTHESIA DISPOSABLE LARYNGOSCOPE

LED cold light source, strong light, matching handle

Cat.No.	Specification
LB371160	160# Curved
LB371140	140# Curved
LB371120	120# Curved
LB371095	95# Curved
LB371130	130# Straight
LB371100	100# Straight
LB371075	75# Straight
LB371000	Handle

LB3720 DISPOSABLE LARYNGOSCOPE SET

ABS material anesthesia laryngoscope
LED cold light source, strong light

Cat.No.	Components
LB3720	storage box
	curved 95mm 1pc
	curved 120mm 1pc
	curved 140mm 1pc
	curved 160mm 1pc
	straight 75mm 1pc
	straight 100mm 1pc
	straight 130mm 1pc
Handle 1pc	



VENTILATOR CIRCUIT

Disposable anesthesia breathing circuit match use with anesthesia machine or breathing machine, as a use of piping anesthesia gases, oxygen and other medical gases into patient. This product is made by non-toxic and smell-less material PP and PE, with the characteristic of good elasticity, flexibility and press tightness



LB431 DISPOSABLE ANESTHESIA BREATHING CIRCUIT-NORMAL

Cat.No.	Size
LB431015	Adult 1.5m
LB431018	Adult 1.8m
LB431020	Adult 2.0m
LB431115	Pediatric 1.5m
LB431118	Pediatric 1.8m
LB431120	Pediatric 2.0m



LB431E DISPOSABLE ANESTHESIA BREATHING CIRCUIT-EXPANDABLE

Cat.No.	Size
LB431015E	Adult 1.5m
LB431018E	Adult 1.8m
LB431020E	Adult 2.0m
LB431115E	Pediatric 1.5m
LB431118E	Pediatric 1.8m
LB431120E	Pediatric 2.0m


LB431R DISPOSABLE REINFORCED ANESTHESIA BREATHING CIRCUIT

Cat.No.	Size
LB431015R	Adult 1.5m
LB431018R	Adult 1.8m
LB431115R	Pediatric 1.5m
LB431118R	Pediatric 1.8m



LB432 DISPOSABLE ANESTHESIA BREATHING CIRCUIT (WITH WATERTRAP)-NORMAL

Cat.No.	Size
LB432018	Adult 1.8m, no holes, 2 watertraps
LB432118	Pediatric 1.8m, no holes, 2 watertraps



LB432R DISPOSABLE ANESTHESIA BREATHING CIRCUIT (WITH WATERTRAP)-REINFORCED

Cat.No.	Size
LB432018R	Adult 1.8m, no holes, 2 watertraps
LB432118R	Pediatric 1.8m, no holes, 2 watertraps



LB4330 BREATHING CIRCUIT-COAXIAL

Cat.No.	Size
LB433015	Adult 1.5m
LB433018	Adult 1.8m



LB4910
ANESTHESIA RESERVOIR - LATEX FREE

Cat.No.	Size
LB491000	0.5L
LB491001	1L
LB491002	2L
LB491003	3L



LB4920
ANESTHESIA RESERVOIR - LATEX

Cat.No.	Size
LB492000	0.5L
LB492001	1L
LB492002	2L
LB492003	3L



LB4930
ANESTHESIA RESERVOIR - SILICONE

Cat.No.	Size
LB493000	0.5L
LB493001	1L
LB493002	2L
LB493003	3L



LB4511
Disposable expandable circuit +HMEF+
Expandable catheter mount +Latex free bag
+ 2 Way connetor for breathing bag

LB4512
Disposable reinforced circuit +HMEF+
Expandable catheter amount +Latex free bag
(Reinforced catheter amount is available)
+ 2 Way connetor for breathing bag




Made of high-quality imported silicone materials and with the method of platinum catalytic technology.
 Patent technology of integrated molding, high strength, not prone to coming off and being separated, good flexibility.
 Can endure high temperature and high pressure disinfection, reusable.
 Joints are made with injection molding with no gas leakage.
 Applicable to all kinds of imported and domestic breathing machine, anesthesia machine. Length could be customized and free options among glass, Y type joint, L-shaped connector, masks, breathing air.



**LB531
 BREATHING TUBE
 FOR ANESRGESIAAPPATUS**
 Size: Adult

**LB532
 BREATHING CIRCUIT
 FOR VENTILATOR**
 Size: Child



**LB4340
 Breathing Circuit-duo Limbo**

Cat.No.	Size
LB434015	Adult 1.5m
LB434018	Adult 1.8m



LB4410R

LB4410

LB4410E

**LB4410
 Catheter Mount**

Cat.No.	Size
LB4410E	Expandable 20cm
LB4410	Normal 20cm
LB4410R	Reinforced 20cm




**LB421007
 Spirometer Filter**

BV for Spirometer
 Filter Efficiency: $\geq 99.9\%$
 Resistance: $< 0.2\text{KPa}$ (at 30L/min)
 Connector: 28mm, 30mm



**LB411001
BV Filter**

Cat.No.	LB411001
Type	Adult/Child
Filtration efficiency(%)	99.9%
Resistance(Pa) 30L/min	less than 100pa
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)



**LB411002
BV Filter**

Cat.No.	LB411002
Type	Adult/Child
Filtration efficiency(%)	99.9%
Resistance(Pa) 30L/min	less than 100pa
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)



**LB421001
HME Filter**

Cat.No.	LB421001
Filtration efficiency(%)	99.999
Airway Resistance	less than 0.14 KPA under 30L/min flow rate.
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)
Humidity Capacity	Adult:50-900ml
Dead Space Capacity	26ml
Weight	34g
Humidification efficiency at VT 250ml	31.8 mg H2O/l



**LB421002
HME Filter**

Cat.No.	LB421002
Filtration efficiency(%)	99.999
Airway Resistance	less than 0.06 KPA under 30L/min flow rate.
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)
Humidity Capacity	Adult:50-900ml
Dead Space Capacity	10ml
Weight	21g
Humidification efficiency at VT 250ml	5.0mg H2O/l

**LB421003
HME Filter**

Cat.No.	LB421003
Filtration efficiency(%)	99.999
Airway Resistance	less than 0.09 KPA under 30L/min flow rate.
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)
Humidity Capacity	Adult:50-900ml
Dead Space Capacity	20ml
Weight	22g
Humidification efficiency at VT 250ml	6.0mg H2O/l



**LB421004
HME Filter**

Cat.No.	LB421004
Type	Adult
Filtration efficiency(%)	99.999
Airway Resistance	less than 0.3KPA under 30L/min flow rate.
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)
Humidity Capacity	Adult:200-1500ml;
Dead Space Capacity	Adult:60ml
Weight	Adult;29g
Humidification efficiency at VT 250ml	5.0mg H2O/l



**LB421005
HME Filter**

Cat.No.	LB421005
Type	Pediatric
Filtration efficiency(%)	99.999
Airway Resistance	less than 0.3KPA under 30L/min flow rate.
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)
Humidity Capacity	Pedia:90-1500ml
Dead Space Capacity	Pedia:30ml
Weight	Pedia: 22g




**LB421006
HME Filter**

Cat.No.	LB421006
Filtration efficiency(%)	99.999
Airway Resistance	less than 0.3KPA under 30L/min flow rate.
Standard Connector	Fitting 15 Male/15 Female - 22 Male (mm)
Humidity Capacity	Adult:200-1500ml
Dead Space Capacity	Adult:60ml
Weight	Adult;29g



Standard Endotracheal Tube (Oral/Nasal Endotracheal Tube)

Endotracheal tube is a temporarily breathing device that is inserted into patient's trachea in order to maintain an open air passage and prevent aspiration of mucus and gastric contents. It is designed for general anaesthesia , intensive care and emergency case for airway management and mechanical ventilation.



LB5010C High Volume, Low Pressure Cuff

Cat.No.	Size ID(mm)
LB501025C	2.5
LB501030C	3.0
LB501035C	3.5
LB501040C	4.0
LB501045C	4.5
LB501050C	5.0
LB501055C	5.5
LB501060C	6.0

Cat.No.	Size ID(mm)
LB501065C	6.5
LB501070C	7.0
LB501075C	7.5
LB501080C	8.0
LB501085C	8.5
LB501090C	9.0
LB501095C	9.5
LB501100C	10.0



LB5010 Uncuffed

Cat.No.	Size ID(mm)
LB501025	2.5
LB501030	3.0
LB501035	3.5
LB501040	4.0
LB501045	4.5
LB501050	5.0
LB501055	5.5
LB501060	6.0

Cat.No.	Size ID(mm)
LB501065	6.5
LB501070	7.0
LB501075	7.5
LB501080	8.0
LB501085	8.5
LB501090	9.0
LB501095	9.5
LB501100	10.0

Endotracheal Tube Oral Preformed

Oral preformed to reduce tube kinking during oral intubation
Used for nasal surgery, ophthalmic surgery, facial surgery.

LB5020C High Volume Low Pressure Cuff

Cat.No.	Size ID(mm)
LB502030C	3.0
LB502035C	3.5
LB502040C	4.0
LB502045C	4.5
LB502050C	5.0
LB502055C	5.5
LB502060C	6.0
LB502065C	6.5

Cat.No.	Size ID(mm)
LB502070C	7.0
LB502075C	7.5
LB502080C	8.0
LB502085C	8.5
LB502090C	9.0
LB502095C	9.5
LB502100C	10.0



LB5020 Uncuffed

Cat.No.	Size ID(mm)
LB502020	2.0
LB502025	2.5
LB502030	3.0
LB502035	3.5
LB502040	4.0
LB502045	4.5
LB502050	5.0
LB502055	5.5

Cat.No.	Size ID(mm)
LB502065	6.5
LB502070	7.0
LB502075	7.5
LB502080	8.0
LB502085	8.5
LB502090	9.0
LB502095	9.5
LB502100	10.0



Endotracheal Tube Nasal Preformed

Nasal preformed to reduce tube kinking during nasal intubation. Used for oral surgery, maxillofacial surgery.



LB5030C High Volume Low Pressure Cuff

Cat.No.	Size ID(mm)
LB503030C	3.0
LB503035C	3.5
LB503040C	4.0
LB503045C	4.5
LB503050C	5.0
LB503055C	5.5
LB503060C	6.0
LB503065C	6.5

Cat.No.	Size ID(mm)
LB503070C	7.0
LB503075C	7.5
LB503080C	8.0
LB503085C	8.5
LB503090C	9.0
LB503095C	9.5
LB503100C	10.0

LB5030 Uncuffed

Cat.No.	Size ID(mm)
LB503020	2.0
LB503025	2.5
LB503030	3.0
LB503035	3.5
LB503040	4.0
LB503045	4.5
LB503050	5.0
LB503055	5.5
LB503060	6.5

Cat.No.	Size ID(mm)
LB503065	6.5
LB503070	7.0
LB503075	7.5
LB503080	8.0
LB503085	8.5
LB503090	9.0
LB503095	9.5
LB503100	10.0

Reinforced Endotracheal Tube

Designed for neurosurgical procedures, head and neck procedures, etc. with a spiral wire embedded in the endotracheal tube provides effective resistance to kinking.

LB5040C High Volume, Low Pressure Cuff

Cat.No.	Size ID(mm)
LB504030C	3.0
LB504035C	3.5
LB504040C	4.0
LB504045C	4.5
LB504050C	5.0
LB504055C	5.5
LB504060C	6.5


Cat.No.	Size ID(mm)
LB504065C	6.5
LB504070C	7.0
LB504075C	7.5
LB504080C	8.0
LB504085C	8.5
LB504090C	9.0
LB504095C	9.5

LB5040 Uncuffed

Cat.No.	Size ID(mm)
LB504030	3.0
LB504035	3.5
LB504040	4.0
LB504045	4.5
LB504050	5.0
LB504055	5.5
LB504060	6.5

Cat.No.	Size ID(mm)
LB504065	6.5
LB504070	7.0
LB504075	7.5
LB504080	8.0
LB504085	8.5
LB504090	9.0
LB504095	9.5





Tracheostomy Tube

Tracheostomy tube is a hollow tube, with or without a cuff, that is electively inserted directly into the trachea through a surgical incision. Or with a wire-guided progressive dilatation technique in case of emergency. Lightweight tube made of medical grade pvc for patient comfort. Full length radio-opaque line for detection of correct placement.

LB5110C High Volume Low Pressure Cuff

Cat.No.	Size ID(mm)	Cat.No.	Size ID(mm)
LB511030C	3.0	LB511070C	7.0
LB511035C	3.5	LB511075C	7.5
LB511040C	4.0	LB511080C	8.0
LB511045C	4.5	LB511085C	8.5
LB511050C	5.0	LB511090C	9.0
LB511055C	5.5	LB511095C	9.5
LB511060C	6.5	LB511100C	10.0
LB511065C	6.5		

LB5110 Uncuffed

Cat.No.	Size ID(mm)	Cat.No.	Size ID(mm)
LB511030	3.0	LB511065	6.5
LB511035	3.5	LB511070	7.0
LB511040	4.0	LB511075	7.5
LB511045	4.5	LB511080	8.0
LB511050	5.0	LB511085	8.5
LB511055	5.5	LB511090	9.0
LB511060	6.5	LB511095	9.5
		LB511100	10.0



LB5210, LB5220 Endotracheal Tube Introducer (Bougie)

Smooth elbow design of tip can assist intubation, decreasing the damage of soft tissue. Scale mark indicates distance to the tip. Recommend for difficult intubation and replacing an endotracheal tube. Available with solid type and hollow type. Hollow type can facilitate delivery of oxygen in emergency.

Cat.No.	Size ID(Fr)	Type
LB521006	6	Solid type
LB521010	10	Solid type
LB521015	15	Solid type
LB522010	10	Hollow type
LB522015	15	Hollow type



LB5310 Intubating Stylet

Made of malleable aluminum covered with plastic sheath. Soft rounded tip for minimizing patient trauma during intubation. Flexible use for assistance of intubation.

Cat.No.	Size ID(Fr)	Appliable ETT	Length(mm)
LB531006	6	2.5-4.5	302
LB531010	10	4.0-6.0	382
LB531012	12	5.5-6.5	382
LB531014	14	6.0-10.0	382






LB5510 Tracheostomy Tube Holder

Cat.No.	Size
LB551001	Adult :430x25mm
LB551002	Pediatric :250x20mm
LB551003	Adult :450x35mm



Surgical Tube

all the materials we use are medical grade. This product is made by non-toxic and smell-less materials, with the characteristic of good elasticity and flexibility.



**LB6100
Suction Yankauer Handle**

Size	Type
LB6101	Crown tip with vent
LB6102	Crown tip without vent
LB6103	Flat tip with vent
LB6104	Flat tip without vent



LB6201



LB6202



LB6203

**LB6200
Connector**

Supplied with universal female connector, can connect with smooth suction yankauer handle.

Size	Type	Color
LB6201	Long Connector	Blue, Green or Transparent
LB6202	Short Connector	Light Blue or Navy Blue
LB6203	Trumpet Connector	Transparent or Blue

Connecting Tube

Made of transparent material for better visualization
Striated walls of tube provide superior strength
and anti- Kinking

Size	ID(mm)	OD(mm)
1/4"	6.5±0.2	9.7±0.2
3/16"	4.7±0.2	8.2±0.2

LB60

Connecting Tube without Yankauer Handle

Supplied with Universal female connector,
can connect with smooth suction yankauer handle.

Code	Suction Tube (Inch)	Connector	Length
LB6010118	1/4	Long connector	1.8m
LB6010220	1/4	Short connector	2m
LB6010336	1/4	Trumpet connector	3.6m
LB6020118	3/16	Long connector	1.8m
LB6020220	3/16	Short connector	2m
LB6020336	3/16	Trumpet connector	3.6m

Length can be customized



LB61

Connecting Tube with Yankauer Handle

Made from Non-toxic PVC, clear and soft
Large Lumen resists clogging and the transparency
Allows Clear visualization of liquids
Crown tip, with/without hole or plain tip
with/without hole.
Connecting tube size: 1/4"x1.8m
1/4"x3.6m, 3/16"x1.8m, 3/16"x3.6m
Packing: Packed in individual blister
packing or Double packing, 20pcs/carton



Code	Suction Tube (inch)	Yankauer Handle	Hole	Connector	Length
LB6111118	1/4	Crown tip	Y	Long connector	1.8m
LB6111218	1/4	Crown tip	Y	Short connector	1.8m
LB6111318	1/4	Crown tip	Y	Funnel connector	1.8m
LB6112118	1/4	Crown tip	N	Long connector	1.8m
LB6112218	1/4	Crown tip	N	Short connector	1.8m
LB6112318	1/4	Crown tip	N	Funnel connector	1.8m
LB6113118	1/4	Flat tip	Y	Long connector	1.8m
LB6113218	1/4	Flat tip	Y	Short connector	1.8m
LB6113318	1/4	Flat tip	Y	Funnel connector	1.8m
LB6114118	1/4	Flat tip	N	Long connector	1.8m
LB6114218	1/4	Flat tip	N	Short connector	1.8m
LB6114318	1/4	Flat tip	N	Funnel connector	1.8m
LB6121118	3/16	Crown tip	Y	Long connector	1.8m
LB6121218	3/16	Crown tip	Y	Short connector	1.8m
LB6121318	3/16	Crown tip	Y	Funnel connector	1.8m
LB6122118	3/16	Crown tip	N	Long connector	1.8m
LB6122218	3/16	Crown tip	N	Short connector	1.8m
LB6122318	3/16	Crown tip	N	Funnel connector	1.8m
LB6123118	3/16	Flat tip	Y	Long connector	1.8m
LB6123218	3/16	Flat tip	Y	Short connector	1.8m
LB6123318	3/16	Flat tip	Y	Funnel connector	1.8m
LB6124118	3/16	Flat tip	N	Long connector	1.8m
LB6124218	3/16	Flat tip	N	Short connector	1.8m
LB6124318	3/16	Flat tip	N	Funnel connector	1.8m

Length can be customized




Plain Type



Cap-Cone



Y-Type



T-Type



Y-Type Transparent

Suction Tube

Made of medical grade PVC
Color coded for size identification
Open distal end with 2 side holes
Different connector can be choose.
Length: 50cm

Size ID (Fr)	Cat.No.:LB6210 (Plain Type)	Cat.No.:LB6220 (Cap-Cone)	Cat.No.:LB6230 (Y-Type)	Cat.No.:LB6240 (T-Type)	Cat.No.:LB6250 (Y-Type Transparent)	Color
Fr6	LB621006	LB622006	LB623006	LB624006	LB625006	Light Green
Fr8	LB621008	LB622008	LB623008	LB624008	LB625008	Blue
Fr10	LB621010	LB622010	LB623010	LB624010	LB625010	Black
Fr12	LB621012	LB622012	LB623012	LB624012	LB625012	White
Fr14	LB621014	LB622014	LB623014	LB624014	LB625014	Green
Fr16	LB621016	LB622016	LB623016	LB624016	LB625016	Orange
Fr18	LB621018	LB622018	LB623018	LB624018	LB625018	Red
Fr20	LB621020	LB622020	LB623020	LB624020	LB625020	Yellow
Fr22	LB621022	LB622022	LB623022	LB624022	LB625022	Purple

Feeding Tube

Made of medical grade PVC
Color coded for size identification
Closed distal end with 2 side holes
With X-ray or without X-ray is available.
Packing in individual PE bag or Blister pack
LB6320 Size: FR4-FR8: 40cm
LB6310 Size: FR10-20: 120cm


Size ID (Fr)	Cat.No.:LB6310 (Type I)	Cat.No.:LB6320 (Type II)	Color
Fr4	LB631004	LB632004	Red
Fr6	LB631006	LB632006	Light Green
Fr8	LB631008	LB632008	Blue
Fr10	LB631010	LB632010	Black
Fr12	LB631012	LB632012	White
Fr14	LB631014	LB632014	Green
Fr16	LB631016	LB632016	Orange
Fr18	LB631018	LB632018	Red
Fr20	LB631020	LB632020	Yellow



LB6310



LB6320





LB6220S Silicone Suction Catheter

Material: Silicone
Size: Fr6 ~ Fr16, Length: 50cm
Open distal end with 1 side hole
Applicable to sputum suction in various clinical departments

Volume: 25ml
Length: 35cm

Size ID(Fr)	Cat.No.:LB6112 (Plain Type)	Cat.No.:LB6122 (Cap-Cone)
Fr6	LB611206	LB612206
Fr8	LB611208	LB612208
Fr10	LB611210	LB612210
Fr12	LB611212	LB612212
Fr14	LB611214	LB612214





**LB6410
Stomach Tube**

Made of medical grade PVC
X-ray opaque line throughout the length
Open distal end with 4 side holes
Length: 125cm or 120cm
Size : Fr6 ~Fr24



**LB6420
Ryle's Stomach Tube**

Made of medical grade PVC
X-ray opaque line throughout the length
Closed distal end with 4 side holes
Length: 125m
Stainless steel balls are sealed in the tip end;
It is easier to indwell the tube at a certain position.

Size ID(Fr)	Cat.No.	Color
Fr6	LB642006	Light Green
Fr8	LB642008	Blue
Fr10	LB642010	Black
Fr12	LB642012	White
Fr14	LB642014	Green
Fr16	LB642016	Orange
Fr18	LB642018	Red
Fr20	LB642020	Yellow
Fr22	LB642022	Purple
Fr24	LB642024	Dark blue



**LB6430
Levin Nasogastric tube**

Made of medical grade PVC
With cap
X-ray opaque line throughout the length
Open distal end with 4 side holes
Length: 125cm or 120cm
Size : Fr6 ~Fr24

**LB6410S
Silicone Stomach Tube**

Made of 100% medical grade silicone,
the tube is soft and clear
X-ray opaque line through the whole tube
Closed distal end with 4 side holes
Clear depth marks on the tube.

Size ID(Fr)	Cat.No.
Fr6	LB641006S
Fr8	LB641008S
Fr10	LB641010S
Fr12	LB641012S
Fr14	LB641014S
Fr16	LB641016S
Fr18	LB641018S




**LB6420S
Silicone Ryle's Stomach Tube**

Made of 100% medical grade silicone
The tube is soft and clear
X-ray opaque line through the whole tube
Closed distal end with 4 side holes
Clear depth marks on the tube.
Stainless steel balls are sealed in the tip end;
It is easier to indwell the tube at a certain position.

Size ID(Fr)	Cat.No.
Fr10	LB642010S
Fr12	LB642012S
Fr14	LB642014S
Fr16	LB642012S
Fr18	LB642014S
Fr16	LB642016S
Fr18	LB642018S






LB6510 Rectal Tube

Made form non-toxic, non-irritant soft PVC
Color coded for size identification
Closed distal end with 2 side holes
With X-ray or without is available
Standard tube length: 40cm

Size ID(Fr)	Cat.No.	Color
Fr18	LB651018	Red
Fr20	LB651020	Yellow
Fr22	LB651022	Purple
Fr24	LB651024	Dark Blue
Fr26	LB651026	Grey
Fr28	LB651028	Green



LB6610 Nelaton Tube

Made of medical grade PVC
Closed distal end with 2 side holes
With X-ray or Without is available
Length: 40cm for male, and 20cm for female

Size ID(Fr)	Cat.No. :LB6610 (Male)	Cat.No. :LB6620 (Female)	Color
Fr6	LB661006	LB662006	Light Green
Fr8	LB661008	LB662008	Blue
Fr10	LB661010	LB662010	Black
Fr12	LB661012	LB662012	White
Fr14	LB661014	LB662014	Green
Fr16	LB661016	LB662016	Orange
Fr18	LB661018	LB662018	Red
Fr20	LB661020	LB662020	Yellow
Fr22	LB661022	LB662022	Purple
Fr24	LB661024	LB662024	Dark blue

LB6810 Closed Suction Catheter

- Sputum suction pipe sealing film produced by imported resin, strong toughness, feel is exquisite and conducive to health care workers
- Sputum suction tube body design scale mark line, facilitate throughout the sputum suction effect
- Closed suction phlegm tube adopts full closed type design, can be used to sync with breathing machine, at the time of sputum suction oxygen to sustainable, convenient and safe operation, to avoid cross-infection
- Specification Fr6-Fr18
- 24 hours

Size ID(Fr)	Cat.No.	Specification
6	LB681006	Tube diameter 2mm, length 450mm, pediatric
8	LB681008	Tube diameter 2.6mm, length 450mm, pediatric
10	LB681010	Tube diameter 3.3mm, length 600mm, pediatric
12	LB681012	Tube diameter 4mm, length 600mm, Adult
14	LB681014	Tube diameter 4.7mm, length 600mm, Adult
16	LB681016	Tube diameter 5.3mm, length 450mm, Adult



**LB4010
Bubble Humidifier**

Oxygen bubble humidifiers are used to provide humidified Oxygen to patients both in hospital or at home
The filter at the end of the inlet tube produces very small Bubbles of gas,thus increasing the contact surface
Providing maximum humidity picked up by the bubbles.
At the same time, the tiny bubbles produce very low noise
Unlike larger bubbles, helping patient rest.
Bottle is supplied with the connector that enables it to be fitted To a fire tree outlet of an oxygen flow meter.
Safety valve at 4 or 6 PSI
Single patient use.
Capacity: 250ml, 500ml



**LB1720
Peak Flow Meter**

It is used to detect the pulmonary for the asthma patients
Material: ABS
Size: 30/50 x 155mm
Packing: 1pc/box, 200pcs/carton



**LB1810
Three Balls Spirometer**



The three balls spirometer is mainly used for normal breathing restored trainings for the patients who has just done splanchnocoele operation

**LB1710
Peak Flow Meter**

It is used to detect the pulmonary for the asthma patients
Material: AS
Size: Adult/Child
Packing: 1pc/box, 50box/carton

