


# Surgical & Examination Gloves



# Surgical & Examination Gloves

**Grip options:**  
Available in smooth or  
micro-roughened finish.

**Optimal finger dimensions:**  
Designed to reduce bagginess  
at the fingertips and reduce  
excessive tightness.

**Comfortable palm design:**  
Designed to enhance comfort and  
reduce hand fatigue throughout  
even the longest procedures.

**Anatomical finger design:**  
Designed to replicate the natural,  
curved shape of the hand.

**Innovative cuff design:**  
Designed with thicker material to  
reduce cuff tearing and a narrower  
circumference to reduce cuff rolldown.

# Surgical & Examination Gloves



## Sterile Surgical Latex Gloves

Code SLG01

- Beaded Cuff
- Textured Fingers and Palm
- Low Allergenicity
- Powdered or Powder Free
- Size: 6, 6.5, 7, 7.5, 8, 8.5, 9



## Sterile Exam Latex Gloves

Code ELG01

- Beaded Cuff
- Textured Fingers and Palm
- Ambidextrous
- Low Allergenicity
- Powdered or Powder Free
- Size: S, M, L, XL



## Exam Latex Gloves

Code ELG02

- Beaded Cuff
- Ambidextrous
- Smooth Surface
- Low Allergenicity
- Powdered or Powder Free
- Size: S, M, L, XL



## Exam Vinyl Gloves

Code EVG01

- Beaded Cuff
- Clear
- Ambidextrous
- Smooth Surface
- Latex Free
- Powdered or Powder Free
- Size: S, M, L, XL



## Exam Nitrile Gloves

Code ENG01

- Beaded Cuff
- Blue
- Ambidextrous
- Textured Fingers and Palm
- Latex Free
- Powdered or Powder Free
- Size: S, M, L, XL



## Exam Synthetic Gloves

Code ESG01

- Beaded Cuff
- Pretty Stretch
- White or Yellow
- Ambidextrous
- Smooth Surface
- Latex Free
- Powdered or Powder Free
- Size: S, M, L, XL

# Examination Gloves & Plastic Items



## Plastic Gloves

Code PG01

- Smooth or Embossed
- Clear/Blue
- Ambidextrous
- Latex Free
- Size: S, M, L, XL



## Smart Gloves

Code SG01

- High Elasticity, Replacement of Vinyl Gloves
- Improved Temperature Performance in High Heat Resistance
- Clear
- Ambidextrous
- Latex Free
- Size: S, M, L, XL



## PE Gown

Code: PEG01

- Made by Polyethylene material
- Totally Resistance From fluid
- Thumb Loop Style
- Blue Color



## PE Apron

Code: PAP01

- Made by 65%LDPE and 35%LLDPE
- Totally Resistance from Fluid
- Smooth or Embossed
- Clear or Blue Color