




GT028-100
Economic Urinary Drainage Bag

Capacity: 2000ml
Without bottom outlet; with 90cm inlet tube,
with non-return valve,
Sterilization: EO
Packing: 1pc/PE bag, 250pcs/carton
Measurement of carton: 52x32x36cm



GT028-200
Economic Urinary Drainage Bag

Capacity: 2000ml
with bottom outlet; with 90cm inlet tube;
with non-return valve
with pull-push valve
Sterilization: EO
Form of packing: 1pc/PE bag, 250pcs/carton
Measurement of carton: 540x400x270mm



GT028-310
Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom; Latex free
with twist valve; with white sample port;
with anti-reflux valve ; with plastic tie
Inlet tube: O.D 10mm; 100cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 100pcs/carton
Measurement of carton: 540x360x320mm



GT028-320
Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom; Latex free
with T valve; with white sample port;
with anti-reflux valve ; without plastic tie
Inlet tube: O.D 10mm; 100cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 100pcs/carton
Measurement of carton: 540x360x320mm



GT028-210
Economic Urinary Drainage Bag

Capacity: 2000ml
With bottom outlet; with 90cm inlet tube;
with non-return valve, with T valve
Sterilization: EO
Packing: 1pc/PE bag, 250pcs/carton
Measurement of carton: 52x37x32cm



GT028-300
Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: square bottom, latex free, with pull push valve,
with white sample port, with anti-reflux valve, with plastic tie
Inlet tube: O.D 10mm; 100cm length
Sterilization: EO
Packing: 1pc/ PE or blister pack, 100pcs/carton,
Measurement of carton: 54x36x32cm



GT028-330
Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom; Latex free
with T valve; with white sample port; with anti-reflux valve
with plastic tie
Inlet tube: O.D 10mm; 100cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 100pcs/carton
Measurement of carton: 540x360x320mm



GT028-340
Luxury Urinary Drainage Bag

Capacity: 2000ml
Type: Square bottom; Latex free
with T valve; without sample port; with anti-reflux valve
with double ropes
Thickness film: 0.11 mm per ply, Front film is transparent
Back film is white
Inlet tube: O.D 6.5mm; 90cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 100pcs/carton
Measurement of carton: 540x360x320mm



GT028-400

Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: Square bottom; Latex free with side valve; with anti-reflux valve with single plastic hanger and tubing clip
 Inlet tube: O.D 10mm; 110cm length
 Sterilization: EO
 Form of packing: 1pc/Blister pack, 40pcs/carton
 Measurement of carton: 430x320x380mm



GT028-410

Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: Square bottom; Latex free with Single reversal valve; with needle sample port with anti-reflux valve; with double plastic hanger and single rope with blue bed sheet clip
 Inlet tube: O.D 10mm; 100cm length
 Sterilization: EO
 Form of packing: 1pc/Blister pack, 40pcs/carton
 Measurement of carton: 430x320x380mm



GT028-600

Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: pyriform; latex free, with 100% silicone blue outlet and with folding clamp, with needle sample port, with anti-reflux valve, with finger guard clamp
 Inlet tube: O.D 10mm; 100cm length
 Sterilization: EO
 Form of packing: 1pc/PE or blister pack, 20pcs/carton.
 Measurement of carton: 38x22x15cm



GT028-700

Urinary Drainage Bag

Capacity: 2000ml
 Type: Pyriform; Latex free with T valve; with needle sample port; with anti-reflux valve with double plastic hanger and rope; with white bed sheet clamp
 Inlet tube: O.D 10mm; 110cm length
 Sterilization: EO
 Form of packing: 1pc/Blister pack, 40pcs/carton
 Measurement of carton: 430x300x380mm



GT028-500

Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: Square bottom; Latex free with twist valve; with anti-reflux valve; with plastic tie
 Front and back films are transparent
 Inlet tube: O.D 6.5mm; 90cm length
 Sterilization: EO
 Form of packing: 1pc/Blister pack, 100pcs/carton
 Measurement of carton: 450x360x330mm



GT028-510

Luxury Urinary Drainage Bag

Capacity: 2000ml
 Type: semicircular bottom, Latex free with twist valve; with anti-reflux valve; with plastic tie
 Front and back films are transparent
 Inlet tube: O.D 6.5mm; 90cm length
 Sterilization: EO
 Form of packing: 1pc/PE bag, 100pcs/carton
 Measurement of carton: 440x360x240mm



GT028-800

Urinary Drainage Bag

Capacity: 2000ml
 Type: Pyriform; Latex free with T valve; with needle sample port; with anti-reflux valve with blue plastic hanger and rope; with white bed sheet clamp
 Inlet tube: O.D 10mm; 100cm length
 Sterilization: EO
 Form of packing: 1pc/Blister pack, 20pcs/carton
 Measurement of carton: 310x260x220mm



GT028-810

Urinary Drainage Bag

Capacity: 2000ml
 Type: Pyriform; Latex free with T valve; with needle sample port; with anti-reflux valve with white double plastic hanger and rope; with white bed sheet clamp
 Inlet tube: O.D 10mm; 110cm length
 Sterilization: EO
 Form of packing: 1pc/Blister pack, 40pcs/carton
 Measurement of carton: 430x300x380mm



**GT028-900
Luxury Urinary Drainage Bag**

Capacity: 4000ml
with T valve; Latex free;
with anti-reflux drip chamber (one part)
with needle sample port and tubing clip is optional
with anti-reflux valve;
with reinforced double plastic hanger and rope
with blue bed sheet clamp
Inlet tube: O.D 11mm; 110cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 40pcs/carton
Measurement of carton: 430x300x380mm



**GT028-910
Luxury Urinary Drainage Bag**

Capacity: Type A : 2000ml / Type B : 2500ml
with T valve; Latex free;
with anti-reflux drip chamber (one part)
with needle Sample port and tubing clip is optional
with anti-reflux valve;
with reinforced double plastic hanger and rope
with white bed sheet clamp
Inlet tube: O.D 10mm; 110cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 40pcs/carton
Measurement of carton: 430x300x380mm



**GT029-300
Urine Meter Drainage Bag**

Capacity: 2600ml+400ml; Latex free
with T valve on urine bag;
with twist valve on urine meter
with needle sample port; with anti-reflux valve
with reinforced double plastic hanger and rope
with white bed sheet clamp
Inlet tube: O.D 11mm; 130cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 10pcs/carton
Measurement of carton: 480x410x250mm



**GT029-310
Urine Meter Drainage Bag**

Capacity: 2600ml +500ml, Latex free
with T valve on urine bag ; with twist valve on urine meter
with needle Sample port; with anti-reflux valve
with reinforced double plastic hanger and rope
with white bed sheet clamp
Inlet tube: O.D 11mm; 130cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 10pcs/carton
Measurement of carton: 480x410x250mm



**GT028-920
Luxury Urinary Drainage Bag**

Capacity: 2000ml
with T valve; Latex free;
with anti-reflux drip chamber (Three part)
with needle sample port and tubing clip is optional
with anti-reflux valve;
with reinforced double plastic hanger and rope
with blue bed sheet clamp
Inlet tube: O.D 11mm; 110cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 40pcs/carton
Measurement of carton: 430x300x380mm



**GT029-100
Luxury Urinary Drainage Bag**

Capacity: 2000ml
with T valve; Latex free;
with anti-reflux drip chamber (Three part)
with needle sample port and tubing clip is optional
with anti-reflux valve;
with reinforced double plastic hanger and rope
with blue bed sheet clamp
Inlet tube: O.D 11mm; 110cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 40pcs/carton
Measurement of carton: 430x300x380mm



**GT029-320
Urine Meter Drainage Bag**

Capacity: 2000ml +200ml; Latex free
with T valve on urine bag; with twist valve on urine meter
with needle sample port; with anti-reflux valve
with reinforced double plastic hanger and rope
with white bed sheet clamp
Inlet tube: O.D 11mm; 130cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 10pcs/carton
Measurement of carton: 480x410x250mm



**GT029-330
Urine Meter Drainage Bag (Double Chambers)**

Easy to read molded plastic meters are calibrated up to 400 ml urine meter
with T valve on urine bag ; with twist valve on urine meter
with needle sample port; with anti-reflux valve
with reinforced double plastic hanger and rope
with white bed sheet clamp
Inlet tube: O.D 11mm; 130cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 10pcs/carton
Measurement of carton: 480x410x250mm



**GT029-340
Urine Meter Drainage Bag**

Easy to read molded plastic meters are calibrated up to 500 ml urine meter
with 2000ml Luxury Urinary Drainage bag
with anti-reflux drip chamber (one part)
with T valve on urine bag; with needle sample port
with anti-reflux valve;
with four belts & two belts to be choose
Inlet tube: O.D 11mm; 130cm length
Sterilization: EO
Form of packing: 1pc/Blister pack, 10pcs/carton
Measurement of carton: 580x290x235mm



**GT029-400
Leg Bag**

Capacity: 600ml
Type: Square
with blue pull-push valve & blue tip cover
with anti reflux valve; with double high elastic belts
Sterilization: EO
Form of packing: 1pc/PE bag, 48pcs/carton
Measurement of carton:
300x260x150mm With 120cm tube



**GT029-430
Leg Bag**

Capacity: 500 ml
Type: ellipse shape
with T valve & short outlet tube;
with double high elastic belts
Sterilization: EO
Form of packing: 1pc/Blister pack, 100pcs/carton
Measurement of carton: 380x260x230mm



**GT029-410
Leg Bag**

Made of clear ,medical grade PVC
Type:ellipse shape
with green twist valve&green tip cover
with anti reflux valve;
with double high elastic belts
Tube length:30cm,60cm,90cm (length can be customized)
On the back of the bag body can be customized flocking or non-woven fabric
Sterilization:EO
Size available:500ml\600ml\ 750ml\ 1000ml
Packing :Blister bag Or PE Bag



**GT029-410
Leg Bag**

Capacity: 500ml & 1000ml
Type: ellipse shape
with green twist valve & green tip cover
with anti reflux valve; with double high elastic belts
Sterilization: EO
Form of packing: 1pc/PE bag, 48pcs/carton
Measurement of carton: 300x260x150mm



**GT029-420
Leg Bag**

Capacity: 500 ml
Type: ellipse shape
with T valve & short outlet tube;
with anti reflux valve
Sterilization: EO
Form of packing: 1pc/Blister pack, 100pcs/carton
Measurement of carton: 550x400x325mm



**GT029-420
Leg Bag(3 chambers)**

Made of clear ,medical grade PVC
Type:ellipse shape
with green twist valve&green tip cover
with anti reflux valve;
with double high elastic belts
Tube length:30cm,60cm,90cm (length can be customized)
On the back of the bag body can be customized flocking or non-woven fabric
Sterilization:EO
Size available:500ml,750ml
Packing :Blister bag Or PE Bag



**GT029-500
Disposable Pediatric Urine Collector**

1. Medical grade PVC material
2. For single use only. Discard after use
3. Do not store in direct sunlight, at extreme temperatures or in high humidity.
4. Do not use if package is open or damaged
5. Before use, the bottom should be pushed to the top of tight to avoid leakage
Size: General Type (Oval Sharp) 100ml/200ml
For Girl Type (Mickey Mouse sharp) 100ml/200ml
Anti-Reflux type 100ml/200ml
Packing :PE Bag