



GT048-100 Surgical Gloves

Description :sterilized Surgical Gloves
(Sterile By Gamma Ray).
Powdered or powder-free
(material): made Of Natural Latex.
Found As Follows:
Products Technical Datas:
Quality:

Specification (mm):

Size	6mm	6.5mm	7mm	7.5mm	8mm	8.5mm	9mm
Length (min,mm)	260	260	270	270	270	280	280
Width (min,mm)	77±5	83±5	89±5	95±5	105±6	108±6	114±5

Before Aging After Aging
(70'cx168hours)tensile Strength (n) ≥ 12 9
Packing:packed In Cartons And In Good Condition.
Packing: 50Pairs/box 500Pairs/ctn, 67 X 26 X 40cm 15/13 Kgs



GT050-100 Vinyl Gloves

Material: Vinyl
Product Name: Vinyl Gloves
Grade: AQL4.0
General Specificatio

Size	Glove Length	Palm Width	Thickness	Content of Powderful
X-Large	240±10mm	14±3mm	0.08±0.01mm	0.1-0.2g
Large	240±10mm	106±3mm	0.08±0.01mm	0.1-0.2g
Middle	240±10mm	106±3mm	0.08±0.01mm	0.1-0.2g
Small	240±10mm	106±3mm	0.08±0.01mm	0.1-0.2g

We can make different weight as following:

S: 3.0g, 3.5g, 4.0g, 4.5g Etc
M: 3.5g, 4.0g, 4.5g, 5.0g Etc
L: 4.0g, 4.5g, 5.0g, 5.5g Etc
XL: 4.5g, 5.0g, 5.5g, 6.0gEtc
Packing: 100pcs/box, 1000Pcs/ctn, 31.5 X 24.5 X 25cm

GT051-100 Nitrile Examination Gloves

Description :nitril Examination Gloves.
Size: XI Arge Large, Medium, Small
Description: Powdered Or Powder-free
Products Technical Datas:
1) specification (mm):

Size	XL	L	M	S
Length (min,mm)	230±10	230±10	230±10	230±10
Width (min,mm)	110±10	100±10	95±10	80±10
Thickness Single(mm)	≥ 0.08	≥ 0.08	≥ 0.08	≥ 0.08

Physical Properties (mm): 26pcs Sampled.
Before Aging After Aging
(100'cx22hours) 14.0 13.0
Tensile Strength (mpa)
Ultimate Elongation(%) 520 410
Packing: 100pcs/box, 2000pcs/carton,
53x27x36cm 18/16kgs



GT049-100 Latex Examination Gloves

Size: Large, Medium, Small
Description: Powdered Or Powder-free
Material: Made Of Natural Latex
Products Technical Datas:
Specification (mm):

Size	L	M	S
Length (min,mm)	230±10	230±10	230±10
Width (min,mm)	110±10	95±10	80±10
Thickness Single(mm)	≥ 0.08	≥ 0.08	≥ 0.08

Physical Properties (mm): 26pcs Sampled.
Before Aging After Aging
(100'cx22hours)tensile Strength (mpa) 21.0 16.0
Ultimate Elongation(%) 700 500
Packing: 100pcs/box,2000pcs/53x27x32cm 18/16kgs

GT049-101 Gynecology Examination Gloves

Material: Natural Latex
Size: S,M,L
Type: Sterile, Powdered, Powder-free
Packing: Paper-Paper Bag

Dimensions	Standards	
	Ideal Medical Glove	En455
Length(mm)	450mm	400mm
Width(mm)		
Small	83±6mm	83±6mm
Medium	95±6mm	95±6mm
Large	108±6mm	108±6mm
Thickness-Single Wall(mm)		
Fingers	>0.2mm	>0.1mm
Palm	>0.16mm	>0.1mm
Cuff	>0.12mm	>0.1mm

Physical Properties:

Property	Standards	
	Ideal Medical Gloves	ASTM D6319(00)
Elongation at break(%)		
Before Aging	Min 750%	Min 700%
After Aging	Min 550%	Min 500%
Tensile Strength(MPa)		
Before Aging	Min 25MPa	Min 24MPa
After Aging	Min 19MPa	Min 18MPa

GT052-100 PE Gloves

MATERIAL: Recycled or virgin polyethylene.
Thickness: 0.018mm, 0.02mm, 0.022mm etc.
Size: L, M, S
Type: HDPE or LDPE, plain or embossed
Color: white. Blue, pink etc.
Different thickness, sizes,
color and packings is available.

GT053-100 Latex Powder-free Finger Cots And Powdered

Latex Powder-free Anti-static Finger Cots.
Designed With A Proprietary Anti-static Agent That Yields Average Surface Resistance Of 108-11Ω Ohms Per Square Surface Unit.
It Virtually Eliminate Electro Static Discharge.
These Finger Cots Are Clean Processed And To Protect Particulates, Chemicals And Solvents With Skin.
Type: Smooth Surface With Rolled, Unrolled Or Cut Is Available.
Specification Size ± 3 mm

	L	M	S
Length(mm)	65	65	65
Width(mm)	29	27	25
Thickness(mm)	0.13	0.13	0.13

Packing: 100Pcs/bag, 30000Pcs/ ctn, 45 X 45 X 30cm



GT060-101 Two Fingers Cots

1.2gm/pc, white
Package: 1pc/bag, 100bag/inner bag,
5000/carton
Carton size: 400x320x27mm

